

PT. SHAMROCK MANUFACTURING CORPORATION

Manufacturer of Latex & Nitrile Gloves

Ji, Rava Medan - Namorambe PS, IV Km. 9 Kab. Delí Serdang - Sumut - Indonesia Tel: (62-61) 7030008; Fax: (62-61) 7030007

K012786

Page Numbers 1 of 2

(1) Name of applicant

: RUDY SALIM

Address

SHAMROCK Manufacturing Company

Jl. Raya Medan - Namorambe PS IV Kabupaten Deli Serdang - Indonesia

Fax No.

Phone No. : 62-61-7030008 : 62-61-7030007

Contact person in U.S.A.

: Emmy Tjoeng

Fax No.

: 626-913-1498

(2) Device details

Trade Name

: Powder free Nitrile Neoprene Examination Gloves

Classification Name

: Powder free Nitrile Neoprene Examination Gloves

(3) Product Code

: 80 LZA

(4) Equivalent device legally

marketed

: Class I Examination Gloves 80 LZA

meeting ASTM D6319-00a

(5) Intended use

: Powder free Nitrile Neoprene Examination Glove is a disposable device intended for medical purpose that is worn on examiner's hand to prevent contamination between patient and examiner.

(6) Technological characteristic of the gloves.

1 2.	Dimensions				
	Sizes	Small	Medium	Large	X-Large
	Length mm (min.)	220	230	230	230 .
	Palm Width mm Thickness	80±10	95±10	111±10	120±10
	1. Cuff mm (min)	0.05	0.05	0.05	0.05
	2. Palm mm(min)	0.05	0.05	0.05	0.05
	3. Finger Tip mm	0.05	0.05	0.05	0.05
b.	Physical Properties				•
	·,		Before ageing		After agoing at 70°C 168 hrs.
	Tensile Strength Ultimate Elongation		: 14 MPa (min) : 500 % (min.)		14 MPa (min) 400 % (min.)

- (7) Performance data is the same as mentioned immediately above.
- (8) Clinical date is not needed for gloves or for most devices cleared by the 510 (K) process.
- (9) Non-clinical data We certify that the gloves meet or exceed ASTM D6319-00a Standard. Meets FDA pinhole requirement. Meets labeling claim.



OCT - 1 2001

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Shamrock Manufacturing Company C/O Ms. Emmy Tjoeng 889 Sotuh Azusa Avenue City Of Industry, California 91748

Re: K012786

Trade/Device Name: Powder Free Nitrile Neoprene Examination Gloves (Purple)

Regulation Number: 880.6250

Regulation Name: Patient Examination Gloves

Regulatory Class: I Product Code: LZA Dated: August 16, 2001 Received: August 20, 2001

Dear Ms, Tjoeng

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration

and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4618. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sinderely yours

Timothy A. Ulatowski

Director

Division of Dental, Infection Control and General Hospital Devices Office of Device Evaluation Center for Devices and Radiological Health

GMP.

PT. SHAMROCK MANUFACTURING CORPORATION

Manufacturer of Latex & Nitrile Gloves

Ji. Raya Medan - Namorambe PS. IV Km. 9 Kab. Deli Serdang - Sumut - Indonesia

Tel: (62-61) 7030008 ; Fax : (62-61) 7030007

K012786

ANNEXURE II

INDICATION FOR USE

Applicant

: PT. SHAMROCK MANUFACTURING CORPORATION

Device Name

: Powderfree Nitrile Neoprene Examination Gloves

Indication for use

Powderfree Nitrile Neoprene Examination Glove is a disposable device intended for medical purpose that is worn on examiner's hand to prevent contamination between patient and examiner.

(signature)

RUDY SALIM

(Type Name)

(date)

(Division Sign-Off)

Division of Dental, Infection Control,

and General Hospital Devices

510(k) Number _