

K013627

NOV 16 2001

**510(K) SUMMARY OF SAFETY AND EFFECTIVENESS**

This summary of safety and effectiveness is provided as part of this Premarket Notification in compliance with 21 CFR, Part 807, Subpart E, Section 807.92.

**1) Submitter's name, address, telephone number, contact person:**

Bob DePalma (714) 889-3070  
Regulatory Affairs  
Medison America, Inc.  
11075 Knott Ave.  
Cypress, CA 90630  
Telephone : 714 - 889 - 3000  
Facsimile : 714 - 889 - 3079  
Email : bdepalma@medison.com

Prepared August 22, 2001

**2) Name of the device, including the trade or proprietary name if applicable, the common or usual name, and the classification name, if known:**Common/Usual Name:

Diagnostic Ultrasound System and Accessories

Proprietary Name:

SA8000 Diagnostic Ultrasound System and Transducers.

<u>Classification Names:</u>	<u>FR Number</u>	<u>Product Code</u>
Ultrasound Pulsed Echo Imaging System	892.1560	90-IYO
Ultrasonic Pulsed Doppler Imaging System	892.1550	90-IYN
Diagnostic Ultrasound Transducer	892.1570	90-ITX

3) **Identification of the predicate or legally marketed device:**

Medison America, Inc. believes that SA 8000 ultrasound system is substantially equivalent to the currently marketed SA 9900 ultrasound system (K002185) and SA8800/HDI1500 ultrasound system (K974269)

4) **Device Description:**

The SA8000 scanner is a general purpose, mobile, software controlled, diagnostic ultrasound system. Its function is to acquire ultrasound data and to display the data as B-mode, M-mode, Color-Flow Doppler, Continuous (CW) Doppler, Pulsed (PW) Doppler and Power Doppler, or as a combination of these modes. The SA8000 also gives the operator the ability to measure anatomical structures and offers analysis packages that provide information that is used to make a diagnosis by competent health care professionals. The SA8000 has real time acoustic output display with two basic indices, a mechanical index and a thermal index, which are both automatically displayed.

Seventeen different models of transducers are available and any three may be connected at the same time. In addition to the initial operational settings for each transducer preprogrammed in the system, user-customized parameter settings for each transducer may be inserted by the operator and stored for recall as needed via the system control panel. Customization includes transmit focusing, filtering, image enhancement processing, dynamic window curve selection. Controls are also provided to select display format (single and various combinations), to activate zoom features, and to utilize the cine loop function. More detailed explanations of these functions and controls are included in Chapter 2 of the Operator manual, and in the software/firmware documentation included in this 510(k) Notification.

The SA8000 uses digital beamforming technology, and supports a variety of Linear, Convex, Phased Array and Static probes for a wide variety of applications. It is an ultrasound scanner, which provides high resolution, high penetration performance, and various measurement functions. Probes are supported in frequencies from 1.0 MHz to 20.0 MHz. These probes can be applied to a variety of clinical applications such as fetal, abdominal, pediatric, small organ, neonatal cephalic, adult cephalic, cardiac, trans-rectal, trans-vaginal, peripheral-vascular, muscular-skeletal, and Intra-operative. The same clinical uses were cleared for the predicate devices, SA9900 (K002185) and SA8800/HDI1500 (K974269)

The system can be used to measure distances and calculate areas, circumferences and volumes, as well as calculate the expected date of delivery by using BPD (biparietal diameter), OFD (occipitofrontal diameter), HC (head circumference), AC (abdominal circumference), AD (abdominal diameter), FL (femur length), CRL (crown rump length), APTD (anteroposterior trunk diameter), TTD (transverse trunk diameter), GS (gestational sac), LMP (last menstrual period.), Cardiac Analysis (volume by area/length, Simpson biplane and single plane, M-mode analysis, Doppler: peak and mean gradients, pressure half time, E/A ratio and continuity equation) and Vascular Analysis (resistive index, pulsatility index, % stenosis, ICA/CCA ratio, Volume flow).

Biopsy guidelines are provided on screen to assist in the collection of tissue samples, using biopsy guide adapters offered as an optional accessory. M-mode uses the scroll display method which has its images flow from the right to the left on the monitor. The SA8000 supports the Cine function (capable of storing up to 256 sequential images) and real-time zoom function to the region-of-interest. The system provides the ability to perform remote viewing of images, without compression, via a DICOM 3.0 compatible output. Management of patient history is possible by image-filing function. High-resolution images are provided by utilizing a technology called digital dynamic receive focusing.

The SA8000 has been designed to meet the following electromechanical safety standards:

- EN 60601-1 (IEC 601-1,) European Norm, Medical Electrical Equipment
- UL 2601-1, Underwriters Laboratories Standards, Medical Electrical Equipment
- C22.2 No. 601.1, Canadian Standards Association, Medical Electrical Equipment
- CEI/IEC 1157:1992, International Electrotechnical Commission, Requirements for the declaration of the acoustic output of medical diagnostic ultrasonic equipment
- EN 60601-1-2 (IEC 60601-1-2,) European Norm, Collateral Standard: Electromagnetic Compatibility
- Compliant with the European Medical Device Directive Certificate issued by TUV.

**5) Intended Use:**

SA8000 intended uses as defined FDA guidance documents are:

- Fetal (includes infertility monitoring of follicle development)
- Abdominal
- Pediatric
- Small Organ
- Neonatal Cephalic
- Adult Cephalic
- Cardiac (Adult, Pediatric)
- Trans-Rectal
- Trans-Vaginal
- Peripheral-Vascular
- Muscular-Skeletal (conventional, superficial)
- Intra-Operative (Neurological, Abdominal, Peripheral vascular)

Typical examinations performed using the system are:

- General abdominal and pelvic studies including organ surveys, assessment, and retroperitoneal cavity studies.
- Study of small parts including breasts, shoulders, thyroid, and the abdominal wall.
- Pediatric scans of organs and bony structures.
- Peripheral vascular applications including carotid arteries, legs, arms, feet, and penile artery.
- Monitoring procedures for infertility studies (other than in vitro fertilization).
- First, second and third trimester pregnancy studies.
- Prostate, prostate biopsy guidance, and rectal wall studies.
- Neonatal head studies.
- Trans-cranial studies of middle cerebral arteries, internal carotid artery, and vertebral arteries.
- Cardiac studies in adults and children.
- Biopsy guidance for tissue or fluid sampling.
- Conventional podiatry scans.

**6) Technological Characteristics:**

This device operates identical to the predicate devices in that piezoelectric material in the transducer is used as an ultrasound source to transmit sound waves into the body. Sound waves are reflected back to the transducer and converted to electrical signals that are processed and displayed as 2D and M-mode, Spectral Doppler, Color Doppler, Power Doppler, or 3D images. Transducer patient contact materials are biocompatible.

The device's acoustic output limits are:

All Applications:

	(Maximum Range)
ISPTA	720 mW/cm <sup>2</sup>
MI	1.9

The limits are the same as predicate Track 3 devices.



**DEPARTMENT OF HEALTH & HUMAN SERVICES**

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

NOV 1 6 2001

Medison America, Inc.  
% Mr. Mark Job  
Program Manager  
TÜV Product Service Inc.  
1775 Old Highway 8 NW, Suite 104  
NEW BRIGHTON MN 55112-1891

Re: K013627

Trade Name: SonoAce SA-8000 Diagnostic Ultrasound System  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulation Number: 21 CFR 892.1560  
Regulation Name: Ultrasonic pulsed echo imaging system  
Regulation Number: 21 CFR 892.1570  
Regulation Name: Diagnostic ultrasonic transducer  
Regulatory Class: II  
Product Code: 90 IYN  
Product Code: 90 IYO  
Product Code: 90 ITX  
Dated: November 1, 2001  
Received: November 5, 2001

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the SonoAce SA-8000 Diagnostic Ultrasound System, as described in your premarket notification:

Transducer Model Number

L5-9EC/5-9 MHz/7.5 MHz/40mm Linear Array  
L5-9ER/5-9 MHz/7.5 MHz/50mm Linear Array  
HL5-9ED/5-9 MHz/7.5 MHz/40mm Linear Array  
L5-10ED/5-10 MHz/7.5 MHz/40mm Linear Array  
LI5-9EV/5-9 MHz/7.5 MHz 40mm Linear Array  
C3-7ED/3-7 MHz/4.5 MHz/50R Curved Array  
EC4-9ES/4-9 MHz/6.5 MHz/ 10R Curved Array  
CL4-8EV/4-8 MHz/6.5 MHz/40R Curved Array  
P2-4AM/2-4 MHz/2.5 MHz Phased Array  
P2-5AC/2-5 MHz/3.5 MHz Phased Array  
P3-7AM/3-7 MHz/5.0 MHz Phased Array  
2.0CW/2 MHz/Static CW  
4.0CW/4 MHz/Static CW  
S-VAW3-5/3-5 MHz/3.5 MHz Volume Curved Array  
S-VAW4-7/4-7 MHz/4.5 MHz Volume Curved Array  
S-VDW5-8(B)/5-8 MHz/6.5 MHz/ Volume Curved Array  
S-VNA5-8(B)/5-8 MHz/6.5 MHz Volume Curved Array

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801, please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

If you have any questions regarding the content of this letter, please contact Rodrigo C. Perez at (301) 594-1212.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure(s)



NOV 16 2001

### Section 4.3 INDICATIONS FOR USE

K013627

#### DIAGNOSTIC ULTRASOUND INDICATIONS STATEMENT

510(k) No.:

System: **SA8000 Ultrasound System**Intended Use: **Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:**

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)	P	P	P		P	Note 1	Notes 2, 7, 8
	Abdominal	P	P	P	P1	P	Note 1	Notes 2, 7, 8
	Intra-operative (See Note 6)	P1	P1	P1		P1	Note 1	Notes 8
	Intra-operative (Neuro.)	P1	P1	P1	P1	P1	Note 1	Notes 8
	Laparoscopic							
	Pediatric	P	P	P	P	P	Note 1	Notes 2, 7, 8
	Small Organ (See Note 5)	P	P	P		P	Note 1	Notes 2, 8
	Neonatal Cephalic	P	P	P	P	P	Note 1	Notes 7, 8
	Adult Cephalic	P	P	P	P	P	Note 1	Notes 8
	Trans-rectal	P	P	P		P	Note 1	Notes 2, 8
	Trans-vaginal	P	P	P		P	Note 1	Notes 2, 8
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)	P	P	P		P	Note 1	Notes 2, 8
	Musculo-skel. (Superfic.)	P2	P2	P2		P2	Note 1	Notes 2, 8
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult	P	P	P	P	P	Note 1	Notes 4, 7, 8
	Cardiac Pediatric	P	P	P	P	P	Note 1	Notes 4, 7, 8
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	P	P	P	N	P	Note 1	Note 8
	Other (spec.)				N			Note 5

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

## Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

Nancy C. Brogdon

K013627

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: L5-9EC / 5-9MHz / 7.5 MHz / 40mm Linear Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric	N	N	N		N	Note 1	Note 2, 8
	Small Organ (See Note 5)	N	N	N		N	Note 1	Note 2, 8
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)	N	N	N		N	Note 1	Note 2, 8
	Musculo-skel. (Superfic.)	N	N	N		N	Note 1	Note 2, 8
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	Note 1	Note 8
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269;  
P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K013627

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: L5-9ER / 5-9MHz / 7.5 MHz / 50mm Linear Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric	N	N	N		N	Note 1	Note 2, 8
	Small Organ (See Note 5)	N	N	N		N	Note 1	Note 2, 8
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)	N	N	N		N	Note 1	Note 2, 8
	Musculo-skel. (Superfic.)	N	N	N		N	Note 1	Note 2, 8
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	Note 1	Note 8
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K013627

## DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT

510(k) No.:

System: SA8000 Ultrasound System

Transducer: HL5-9ED / 5-9MHz / 7.5 MHz / 40mm Linear Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric	N	N	N		N	Note 1	Note 2, 8
	Small Organ (See Note 5)	N	N	N		N	Note 1	Note 2, 8
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)	N	N	N		N	Note 1	Note 2, 8
	Musculo-skel. (Superfic.)	N	N	N		N	Note 1	Note 2, 8
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	Note 1	Note 8
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 901.109)



(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number 5013627

Basic Information

Section 4.3, Page 4

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: L5-10ED / 5-10 / 7.5 MHz / 40mm Linear Array

Intended Use: Diagnostic <sup>7.5 MHz</sup>ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric	N	N	N		N	Note 1	Note 2, 8
	Small Organ (See Note 5)	N	N	N		N	Note 1	Note 2, 8
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)	N	N	N		N	Note 1	Note 2, 8
	Musculo-skel. (Superfic.)	N	N	N		N	Note 1	Note 2, 8
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel	N	N	N		N	Note 1	Note 8
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K013627

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: LI5-9EV / 5-9MHz / 7.5 MHz / 40mm Linear Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)	N	N	N		N	Note 1	Notes 8
	Intra-operative (Neuro.)	N	N	N	N	N	Note 1	Notes 8
	Laparoscopic							
	Pediatric							
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)


Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)



Basic Information

(Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number K013627

Section 4.3, Page 6

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: C3-7ED / 3-7MHz / 4.5 MHz / 50R Curved Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)	N	N	N		N	Note 1	Note 2, 7, 8
	Abdominal	N	N	N		N	Note 1	Note 2, 7, 8
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric	N	N	N		N	Note 1	Note 2, 7, 8
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E  
Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Basic Information

*Nancy C. Brogdon*  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K013627

Section 4.3, Page 7

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: EC4-9ES / 4-9MHz / 6.5 MHz / 10R Curved Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric							
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal	N	N	N		N	Note 1	Note 2, 8
	Trans-vaginal	N	N	N		N	Note 1	Note 2, 8
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269;  
P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)


Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)



Basic Information

(Division Sign-Off)

Section 4.3, Page 8

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K013627



**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: CL4-8EV / 4-8MHz / 6.5 MHz / 40R Curved Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)	N	N	N		N	Note 1	Notes 8
	Intra-operative (Neuro.)	N	N	N	N	N	Note 1	Notes 8
	Laparoscopic							
	Pediatric							
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)  
Prescription Use (Per 21 CFR 801.109)

## DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT

510(k) No.:

System: SA8000 Ultrasound System

Transducer: P2-4AM / 2-4 MHz / 2.5 MHz Phased Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal	N	N	N	N	N	Note 1	Note 7
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric							
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic	N	N	N	N	N	Note 1	Notes 7
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult	N	N	N	N	N	Note 1	Note 4, 7
	Cardiac Pediatric	N	N	N	N	N	Note 1	Note 4, 7
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K013627

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: P2-5AC / 2-5 MHz / 3.5 MHz Phased Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal	N	N	N	N	N	Note 1	Note 7
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric							
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic	N	N	N	N	N	Note 1	Note 7
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult	N	N	N	N	N	Note 1	Note 4, 7
	Cardiac Pediatric	N	N	N	N	N	Note 1	Note 4, 7
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K013627

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: P3-7AM / 3-7MHz / 5.0 MHz Phased Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal	N	N	N	N	N	Note 1	Note 7
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric	N	N	N	N	N	Note 1	Note 7
	Small Organ (See Note 5)							
	Neonatal Cephalic	N	N	N	N	N	Note 1	Note 7
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric	N	N	N	N	N	Note 1	Note 4, 7
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number

*Nancy C. Brogdon*

K013627

Basic Information

Section 4.3, Page 12

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: 2.0CW / 2MHz / Static CW

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric							
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic				N			
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult				N			
	Cardiac Pediatric				N			
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel				N			
	Other (spec.)				N			Note 5

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

Basic Information

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K013627

Section 4.3, Page 13

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: 4.0CW / 4MHz / Static CW

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric							
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult				N			
	Cardiac Pediatric				N			
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel				N			
	Other (spec.)				N			Note 5

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269;  
P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)



(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

K013627

Basic Information

Section 4.3, Page 14

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: S-VAW3-5 / 3-5MHz / 3.5 MHz Volume Curved Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)	P2	P2	P2		P2	Note 1	Notes 2, 7, 8
	Abdominal	P2	P2	P2		P2	Note 1	Notes 2, 7, 8
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric	P2	P2	P2		P2	Note 1	Notes 2, 7, 8
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269;

P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)


Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)



(Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number K013627

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: S-VAW4-7 / 4-7MHz / 4.5 MHz Volume Curved Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)	P2	P2	P2		P2	Note 1	Notes 2, 7, 8
	Abdominal	P2	P2	P2		P2	Note 1	Notes 2, 7, 8
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric	P2	P2	P2		P2	Note 1	Notes 2, 7, 8
	Small Organ (See Note 5)	N	N	N		N	Note 1	Notes 2, 7, 8
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*Nancy C. Braden*

K013627



**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: S-VDW5-8(B) / 5-8MHz / 6.5 MHz Volume Curved Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric							
	Small Organ (See Note 5)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal	P2	P2	P2		P2	Note 1	Notes 2, 8
	Trans-vaginal	P2	P2	P2		P2	Note 1	Notes 2, 8
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269; P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number

*Nancy C. Brogdon*  
*K013627*

**DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE STATEMENT**

510(k) No.:

System: SA8000 Ultrasound System

Transducer: S-VNA5-8(B) / 5-8MHz / 6.5 MHz Volume Curved Array

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation (*includes simultaneous B-mode)						
General (Track I only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler*	Combined* (Spec.)	Other (Spec.)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal (See Note 3)							
	Abdominal							
	Intra-operative (Abdominal, vascular)							
	Intra-operative (Neuro.)							
	Laparoscopic							
	Pediatric	P3	P3	P3		P3	Note 1	Note 8
	Small Organ (See Note 5)							
	Neonatal Cephalic	P3	P3	P3		P3	Note 1	Note 8
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Cardiac)							
	Musculo-skel. (Convent.)							
	Musculo-skel. (Superfic.)							
	Intra-luminal							
	Other (spec.)							
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Trans-esophageal (Cardiac)							
	Other (spec.)							
Peripheral Vessel	Peripheral vessel							
	Other (spec.)							

N= new indication; P= previously cleared by FDA in K002185 and K974269; P1= previously cleared by FDA in K974269;

P2= previously cleared by FDA in K002185; P3= previously cleared by FDA in K984639; E= added under Appendix E

Additional Comments:

**Color Doppler includes Power (Amplitude) Doppler (P1, P2, P3)**

Note 1: PWD/Color Doppler, PWD/Power Doppler, CWD/Color Doppler, CWD/Power Doppler (P1, P2, P3)

Note 2: Includes imaging for guidance of biopsy (P1, P2)

Note 3: Includes infertility monitoring of follicle development (P1, P2)

Note 4: Color M-mode (P1)

Note 5: For example: thyroid, parathyroid, breast, scrotum and penis in adult, pediatric and neonatal patients (P1, P2)

Note 6: Abdominal organs and peripheral vessel (P1, P2)

Note 7: Tissue Harmonic Imaging (THI) (P1, P2)

Note 8: 3D Imaging (P1, P2, P3)

Concurrence of CDRH, Office of Device Evaluation (ODE)


Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal,

and Radiological Devices

510(k) Number



K013627

Section 4.3, Page 18