

1 510(k) Summary of Safety And Effectiveness for Model LAD-01

1.1 Submitter Information

Norwood Abbey Limited. 63 Wells Road
Chelsea Heights VIC 3196
Australia

Contact Person: Andrew Ellis
Telephone No.: +61-3-9782 7344

1.2 Device Name

Classification Name: Laser surgical instrument for general and plastic surgery
and in dermatology
Proprietary Name: LAD, model LAD-01

1.3 Predicate Devices

Currently marketed Phoresor II Model 900 (K974855) and the LaserLancet LB100 (K980538) were selected as the predicate devices for this submission.

1.4 Description of the Device

The Norwood Abbey Model LAD-01 is a portable, battery powered, single pulse Er:YAG laser. The radiant energy produced by this laser has a wavelength of 2.94 μm , duration of approximately 300 μs , and a beam diameter or spot-size of 3 mm at the treatment site. The radiation delivered by the device is sufficient to remove the stratum corneum of skin exposed to the treatment.

1.5 Indications for Use

The LAD device is indicated for the rapid production of local dermal anesthesia using topical OTC lidocaine 4%.

1.6 Safety & Effectiveness

An IRB, IDE approved study was conducted over two centers and involving 320 subjects. Full needle insertion of a 23G 5/8" needle was used as a standard pain stimulus. After needle insertion subjects were asked to rate the pain felt using an accepted pain evaluation model; pain scores were tested and analyzed as per the approved protocol.

Results revealed that the model LAD-01 can produce a rapid dermal anaesthetic effect using topical OTC lidocaine 4%.

Analysis of clinical observations from 318 subjects taken approximately 48 hours after treatment revealed no adverse events over those of the placebo treatment.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JAN 08 2003

Norwood Abbey Limited
c/o Mr. Robert T. Handren, Jr.
President
Handren Associates, Inc.
5818 Princess Caroline Place
Leesburg, Florida 34748

Re: K021222

Trade/Device Name: LAD Model LAD-01

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser surgical instrument for use in general and
plastic surgery and in dermatology

Regulatory Class: II

Product Code: GEX

Dated: October 8, 2002

Received: October 10, 2002

Dear Mr. Handren:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Mr. Robert T. Handren, Jr.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97) you may obtain. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink that reads "Miriam C. Provost". The signature is written in a cursive style with a large, stylized 'M' and 'P'.

for Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K021222

Device Name: LAD Model LAD-01

Indications for Use:

The LAD is indicated for the ablation of the outer layer of the skin prior to the application of OTC topical 4% lidocaine cream, for local dermal anesthesia.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Optional Format 3-10-98)

Miriam C. Provost
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K021222