

JUL 24 2002

K021322

510(k) Summary

Pursuant to 512(i)(3)(A) of the Food, Drug and Cosmetic Act, Intra-Lock International is required to submit with this Premarket Notification either an "...adequate summary of any information respecting safety and effectiveness or state that information will be made available upon request of any person." Intra-Lock International chooses to submit a summary of the safety and effectiveness information. The summary is as follows:

Trade Name: Intra-Lock System

Sponsor: Intra-Lock International
1200 North Federal Highway
Suite 200
Boca Raton, FL 33432
Registration No.: 3003631996

Device Generic Name: Dental Implants

Classification: According to Section 513 of the Federal Food, Drug, and Cosmetic Act, the device classification is Class III.

Predicate Devices:

Branemark Dental Implant System	K944683
Branemark MkII Dental Implants	K945398
Lifecore Self Tapping Implants	K950624
Lifecore Stage 1 RBM Implants	K003226
ITI Strauman One Part Implants	K994104

Product Description:

The Intra-Lock International Dental Implant System consists of root form dental implants and restorative components which provide the clinician with screw retained, cement retained and removable overdenture type restorative options. This system also includes prosthetic instrumentation, surgical drills, hand instruments, surgical trays and handpiece adapters for use in the surgical procedures. The implants are sterile packaged and include cover screws and a placement instrument. Prosthetic devices are packaged separately to allow the clinician to choose the appropriate means for restoration after osseointegration. Implant bodies are subjected to a surface treatment of blasting with resorbable blast media and then acid passivation to roughen the surface.

Indications for Use:

The Intra-Lock Implant System has been designed to restore partially or fully edentulous patients. The implants have been designed to be used in either the mandible or maxilla and to support removable or fixed prosthesis, from single tooth replacement to full arch reconstruction.

RECEIVED
JUL 18 2002
Page 35

Safety and Performance:

This submission is a Special 510(k): Abbreviated 510(k) as described in FDA's guidance document entitled "The New 510(k) Paradigm - Alternate Approaches to Demonstrating Substantial Equivalence in Premarket Notifications." In support of this 510(k), Intra-Lock International has provided information to demonstrate conformity with FDA's guidance document entitled: Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Abutments Draft Guidance for Industry and FDA.

Conclusion:

Based on the indications for use, technological characteristics, and comparison to predicate devices, the Intra-Lock Dental Implant System has been shown to be safe and effective for its intended use.

RECEIVED
JUL 19 2011



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUL 24 2002

Mr. Jeffery Sakoff
Director of Operations
Intra-Lock International
1200 N. Federal Highway, Suite 200
Boca Raton, Florida 33432

Re: K021322
Trade/Device Name: Intra-Lock Dental Implant System
Regulation Number: 872.3640
Regulation Name: Endosseous Implant
Regulatory Class: III
Product Code: DZE
Dated: April 19, 2002
Received: April 25, 2002

Dear Mr. Sakoff:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Timothy A. Ulatowski
Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K021322

Device Name: Intra-Lock Dental Implant System

Indications for Use:

The Intra-Lock Implant System has been designed to restore partially or fully edentulous patients. The implants have been designed to be used in either the mandible or maxilla and to support removable or fixed prosthesis, from single tooth replacement to full arch reconstruction.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓
(Per 21 CFR 801.109)

OR

Over-the-Counter Use _____

Robert S. Ketz for Dr. Susan Reamer
(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K021322