

AUG 02 2002

K021455

## 510(K) Summary

### 1. Device Trade Name

6F RAPIDO™ Guiding Catheter

8F RAPIDO™ Guiding Catheter

### 2. Device Common Name

Percutaneous Catheter

### 3. Device Description

The RAPIDO Guiding Catheters are available in two French sizes (6F and 8F). The 6F RAPIDO Guiding Catheter has a standard working length of 69 cm. The 8F RAPIDO Guiding Catheter has as standard working length of 47 cm.

### 4. Intended Use

The Guidant Guiding Catheter is intended to access the coronary venous system, and may be used as a dual-catheter assembly. The catheter serves as a conduit for the delivery of contrast medium and devices, including implantable coronary venous leads, introduced into the coronary venous system.

### 5. Technological Characteristics

Comparisons of the 6F/8F RAPIDO Guiding Catheters and predicate devices show that the technological characteristics such as product performance, design, and intended use are substantially equivalent to the currently marketed predicate devices.

### 6. Performance Data

Testing demonstrated that the 6F/8F RAPIDO Guiding Catheters met the acceptance criteria and performed similarly to the predicate device (EASYTRAK Guiding Catheter). No new safety or effectiveness issues were raised during the testing program. The 6F/8F RAPIDO Guiding Catheter may be considered substantially equivalent to the predicate device.

### 7. Conclusion

The Guidant 6F/8F RAPIDO Guiding Catheters are substantially equivalent to the currently marketed EASYTRAK Guiding Catheter (K021284, May 2, 2002) with regards to intended use and design.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

AUG 02 2002

Guidant Corporation  
c/o Ms. Karen S. Alsop  
Principal Regulatory Affairs Associate  
Cardiac Rhythm Management  
4100 Hamline Avenue North  
St. Paul, MN 55112-5798

Re: K021455

Trade Name: RAPIDO™ Guiding Catheters, 6F and 8F  
Regulation Number: 21 CFR.870.1250  
Regulation Name: Percutaneous Catheter  
Regulatory Class: Class II (two)  
Product Code: DQY  
Dated: May 3, 2002  
Received: May 6, 2002

Dear Ms. Alsop:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4646. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Donna-Bea Tillman".

Donna-Bea Tillman, Ph.D.  
Acting Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## INDICATIONS FOR USE STATEMENT

**Applicant:** Guidant Corporation  
4100 Hamline Avenue North  
St. Paul, MN 55112

**510(k) Number (if known):** *To be assigned by FDA*

**Device Names:** RAPIDO™ Guiding Catheter, 6F and 8F

**Intended Use/Indications for Use:**

The Guidant Guiding Catheter is intended to access the coronary venous system, and may be used as a dual-catheter assembly. The catheter serves as a conduit for the delivery of contrast medium and devices, including implantable coronary venous leads, introduced into the coronary venous system.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)


---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-The-Counter \_\_\_\_\_  
(Optional Format 1-1-96)

  
Division of Cardiovascular & Respiratory Devices  
510(k) Number K021433