

JUL 30 2002

### 3 Summary of Safety and Effectiveness

This summary of 510(k) safety and effectiveness information is submitted in accordance with the requirements of the Safe Medical Devices Act of 1990 as implemented in 21 C.F.R. §807.92.

The submitter of this premarket notification is:

Paul Smolenski  
Regulatory Engineer  
Philips Ultrasound  
3000 Minuteman Road, MS 0135  
Andover, MA 01810  
Tel: (978) 659-3380  
Fax: (978) 975-7324

This summary was prepared on June 21<sup>st</sup>, 2002

The proprietary name of the device is the M2424 Diagnostic Ultrasound System with the 21315 ultrasound transducer. These devices are commonly known as a diagnostic ultrasound system and transducer.

These devices are classified as follows:

90IYN	Ultrasonic Pulsed Doppler Imaging System
90IYO	Ultrasonic Pulsed Echo Imaging System
90ITX	Diagnostic Ultrasound Transducer

As stated in 21 CFR, parts 892.1550, 892.1560 and 892.1570, each of these generic types of devices have been classified as Class II.

The M2424 is a diagnostic ultrasound device. It consists of a system console containing the power supply and electronic circuitry required to generate the image, a display screen, and a connection to the separate transducer. It is substantially equivalent to the generic class of ultrasound systems including previous versions of the Philips Medical Systems M2424 (K002470) system and the Voluson 730 (K003525).

The 21315 transducer is substantially equivalent to the generic class of ultrasound transducers including previously cleared Philips Medical Systems transducers (K002470).

The M2424 system and 21315 transducer function in a manner identical to all ultrasound systems and transducers. The system circuitry generates an electronic voltage pulse which is transmitted to the transducer. In the transducer, a piezo electric array converts the electronic pulse into a ultrasonic pressure wave. When coupled to the body, the

pressure wave transmits through body tissues. The differing acoustic properties of the tissues in the body reflect some of the transmitted energy back to the transducer, where it is converted back to electrical signals, processed, and sent back to the system. In the system, advanced signal processing technologies further convert the returned signals into images of the tissues. The Doppler functions of this system process the Doppler shift frequencies from the echoes of moving targets (such as blood), to detect and graphically display the Doppler shifts of these tissues as flow.

The M2424 system with the 21315 transducer is intended for diagnostic ultrasound imaging and fluid flow analysis.

The M2424 system with the 21315 transducer is substantially equivalent in safety and effectiveness to the predicates identified above:

- Both the predicate device and the M2424 are indicated for the diagnostic ultrasonic imaging and fluid flow analysis.
- Both the predicate device and the M2424 have the same gray-scale and Doppler capabilities.
- Both the predicate device and the M2424 use essentially the same technologies for imaging, Doppler functions and signal processing.
- Both the predicate device and the M2424 have acoustic output levels below the applicable FDA limits.
- Both the predicate device and the M2424 are manufactured under equivalent quality systems.
- Both the predicate device and the M2424 are manufactured of materials with equivalent biosafety. The materials have been evaluated and found to be safe for this application.
- Both the predicate device and M2424 are designed and manufactured to the same electrical and physical safety standards.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

**JUL 30 2002**

Philips Ultrasound, Inc.  
c/o Mr. Mark Job  
Responsible Third Party Official  
TUV Product Service  
1775 Old Highway 8 NW, Suite 104  
NEW BRIGHTON MN 55112-1891

Re: K022303

Trade Name: Philips Ultrasound M2424 Diagnostic Ultrasound System  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulation Number: 21 CFR 892.1560  
Regulation Name: Ultrasonic pulsed echo imaging system  
Regulation Number: 21 CFR 892.1570  
Regulation Name: Diagnostic ultrasound transducer  
Regulatory Class: II  
Product Code: 90 IYN, 90 IYO and 90 ITX  
Dated: July 15, 2002  
Received: July 16, 2002

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the Philips Ultrasound M2424 Diagnostic Ultrasound System, as described in your premarket notification:

Transducer Model Number

21315

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:


Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801, please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

If you have any questions regarding the content of this letter, please contact Rodrigo C. Perez at (301) 594-1212.

Sincerely yours,

  
*for*

Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure(s)

K 02230.

## 4.3.2 Indications for Use Summary

## DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE FORM

System: Philips Ultrasound M2424 Diagnostic Ultrasound System

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track I Only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined	Other
Ophthalmic	Ophthalmic	P		P		P	P	
Fetal Imaging & Other	Fetal	P	P	P	P	P	P	2,3
	Abdominal	P	P	P	P	P	P	2,3
	Intra-operative (Specify) <sup>6</sup>	P	P	P	P	P	P	2,3
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	P	P	P	P	P	P	2,3,4,5
	Small Organ (Specify) <sup>7</sup>	P	P	P	P	P	P	2,3
	Neonatal Cephalic	P	P	P	P	P	P	2,3
	Adult Cephalic	P	P	P	P	P	P	2,3
	Trans-rectal	P	P			P	P	
	Trans-vaginal	P	P			P	P	
	Trans-urethral							
	Trans-esoph. (non-Card.)	P	P	P	P	P	P	2,3
	Intra-luminal							
	Other (Specify)							
Cardiac	Cardiac Adult	P	P	P	P	P	P	2,3,4,5
	Cardiac Pediatric	P	P	P	P	P	P	2,3,4,5
	Trans-esoph. (Cardiac)	P	P	P	P	P	P	2,3
	Other (Specify)							
Peripheral Vessel	Peripheral vessel	P	P	P	P	P	P	2,3,4,5
	Other (Specify)							
	Musculo-skel (conventional)							
	Musculo-skel (superficial)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

## Notes:

1. Combined modes are: B+M, B+M+Color, B+PW, B+PW(CVI), B+PW+Color
  2. Harmonics (Tissue and Contrast) - Previously Cleared
  3. Angio - Previously Cleared
  4. Real-time 3D imaging- New
  5. Real-time Bi-plane imaging- New
  6. Intraoperative applications include: Cardiac, Vascular - Previously Cleared
  7. Small Parts applications include: Thyroid, scrotum, prostate, breast - Previously Cleared
- (PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)
- Concurrence of Center for Devices and Radiological Health, Office of Device Evaluation

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number

K022303

# DIAGNOSTIC ULTRASOUND INDICATIONS FOR USE FORM

System: Philips Ultrasound M2424 Diagnostic Ultrasound System

Transducer: 21315 Phased array transducer

Intended Use: Diagnostic ultrasound imaging and fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track I Only)	Specific (Tracks I & III)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	Other (Specify)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal	N	N	N	N	N	N	2,3,4,5
	Abdominal	N	N	N	N	N	N	2,3,4,5
	Intra-operative (Specify) <sup>1</sup>	N	N	N	N	N	N	2,3,4,5
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	N	N	N	N	N	N	2,3,4,5
	Small Organ (Specify)							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Intra-luminal							
	Other (Specify)							
Cardiac	Cardiac Adult	N	N	N	N	N	N	2,3,4,5
	Cardiac Pediatric	N	N	N	N	N	N	2,3,4,5
	Trans-esoph. (Cardiac)							
	Other (Specify)							
Peripheral Vessel	Peripheral vessel							
	Other (Specify)							
	Musculo-skel (conventional)							
	Musculo-skel (superficial)							

N= new indication; P= previously cleared by FDA; E= added under Appendix E

\*Examples may include: A-mode, Amplitude Doppler, 3-D Imaging, Harmonic Imaging, Tissue Motion Doppler, Color Velocity Imaging

Additional Comments: Other Modes include Amplitude Doppler.

Note 1: Intraoperative uses include abdominal and cardiac (including vessels) applications.

Note 2: Other conventional modes includes Angio

Note 3: Real-time 3D imaging - New


Note 4: Real-time Bi-plane imaging - New

Note 5: Harmonics (Tissue and Contrast)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of Center for Devices and Radiological Health, Office of Device Evaluation

Prescription Use (Per 21 CFR 801.109)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number K02230.3