

SEP 19 2002

1022363

510(K) SUMMARY
(as Required by 21 CFR § 807.92)

A. Submitters Information

Submitter's Name: St. Jude Medical, Inc
Cardiac Surgery Division

Address: St. Jude Medical, Inc.
One Lillehei Plaza
St. Paul, MN 55117

Contact Name William McKelvey
Regulatory Affairs Specialist
St. Jude Medical, Inc.
Bus : (651) 415-7029
Fax : (651) 766-3049
Email: wmckelvey@sjm.com

Submission Prepared: July 18, 2002

B. Device Information

Proprietary Name: SJM® Tailor™ annuloplasty band
(Tailor band) Model TAB-(size)

Common or Usual Name: Flexible Annuloplasty band
"C" ring

Classification: Pre-amendment Class II CFR § 870.3800
Cardiovascular Prosthetic Devices,
Annuloplasty Ring

Predicate Device: St. Jude Medical considers The Tailor
annuloplasty band, model TAB to be
substantially equivalent to the Tailor
annuloplasty ring model TARP.

Device Description The Tailor band is a flexible band fabricated
with a medical grade silicone rubber core,
surrounded by polyester fabric.

Intended Use:

The SJM® Tailor annuloplasty band is indicated for use in repair of a mitral or tricuspid heart valve that is diseased or damaged due to acquired or congenital processes.

C. Comparison of Required Technological Characteristics

SJM considers the Tailor band, model TAB to be substantially equivalent in configuration, function and intended use to the Tailor ring, model TARP. The table below is a comparison of the equivalency characteristics between the two devices.

Characteristic	Equivalency
a. Product Labeling	Substantially Equivalent
b. Intended Use	Identical
c. Physical Characteristics	Substantially Equivalent
d. Anatomical Sites	Identical
e. Target Population	Identical
f. Performance Testing	Identical
g. Safety Characteristics	Identical

D. Summary of Non-Clinical Tests

The testing for the Tailor ring model TARP (predicate) is included in the premarket notification (K014161). The following tests have been performed on the Tailor band model TAB to insure substantial equivalence with the predicate.

1. Physical Validation
 - Band end evaluation
 - Gapping
 - Band release
 - Suture Retention
 - Suture pull
2. Microbiological Verification
 - Biocompatibility evaluation
 - Sterility Assurance evaluation
3. Manufacturing Process Validation



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

SEP 19 2002

St. Jude Medical, Inc.
c/o Mr. William McKelvey
Regulatory Affairs Coordinator
One Lillehei Plaza
St. Paul, MN 55117

Re: K022363
Trade Name: SJM® Tailor™ Annuloplasty Band, Model TAB-(size)
Regulation Number: 21 CFR 870.3800
Regulation Name: Annuloplasty Ring
Regulatory Class: Class II (two)
Product Code: KRH
Dated: September 5, 2002
Received: September 6, 2002

Dear Mr. McKelvey:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4586. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Bram D. Zuckerman", is written over the typed name.

Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): _____


Device Name: SJM® Tailor™ annuloplasty band

Indications for Use:

The SJM® Tailor™ annuloplasty band is indicated for use in the repair of a mitral or tricuspid valve that is diseased or damaged due to acquired or congenital processes.

(Please Do Not Write Below This Line-Continue On Another Page If Another Page Is Needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)


Division of Cardiovascular & Respiratory Devices
510(k) Number K022363

Prescription Use X

or

Over-The-Counter Use _____

(Per 21 CFR 801.109)

(Optional Format 1-2-96)