

FEB 27 2003

## **510 (K) SUMMARY**

Date of Summary: 2-17-03

**Product Name:**

STARTOX (4)

**Sponsor:**

Starplex Scientific, Inc.  
50 Steinway Blvd.  
Etobicoke, Ontario, Canada M9W 6Y3

**Correspondent:**

Fran White  
MDC Associates  
163 Cabot Street  
Beverly, MA 01915

**Substantially Equivalent Devices:**

**Product: Rapid Drug Screen – 9 panel**  
**Manufactured by: American BioMedica Corp**  
**K Number: K012159**

**PRODUCT DESCRIPTION:**

A lateral flow immunoassay for the detection of drugs of abuse.

510k Submission  
Starplex, Inc.  
Request for additional data  
K022388

**INTENDED USE:**

“STARTOX” Drugs of Abuse Screening Test is a one-step lateral flow immunoassay intended for the simultaneous detection of multiple drug analytes in urine. “STARTOX” is intended for use in the qualitative detection of drugs of abuse at the following Substance Abuse Mental Health Services Administration (SAMHSA) recommended levels:

Compound	Abbreviation	Level
Methamphetamine ((+/-)methamphetamine HCl)	METH	1000 ng/ml
Opiates (morphine)	OPI	300 ng/ml
Cocaine (benzoylecgonine)	COC	300 ng/ml
Cannabinoids (11-nor- $\Delta^9$ -THC-9-carboxylic-acid)	THC	50 ng/ml

**STARTOX™ Drugs of Abuse Screening Test provide only a preliminary qualitative test result. Use a more specific alternate quantitative analytical method to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Apply clinical and professional judgment to any drug of abuse test result, particularly when preliminary positive results are obtained.**

**PERFORMANCE CHARACTERISTICS:**

STARTOX™ drugs of abuse screening test detects 4 drugs in human urine at the levels indicated.

STARTOX™ is substantial equivalent to Rapid Drug Screen 9 panel manufactured by American BioMedica Corporation 510k number K002447.

Product performance for Opiates, Cocaine and Cannabinoids was compared Medtox Profile II manufactured by Medtox Diagnostics. Performance of Methamphetamine was compared to RDS 9 manufactured by American BioMedica, Inc. Ninety (90) samples were tested against each drug, 50 negative and 40 positive specimens. Drug status of these samples was tested by EMIT and quantified by GC/MS.

The STARTOX™ Drugs of Abuse Screening Test was compared to a FDA substantially equivalent approved device. STARTOX™ performed substantially equivalent to the legally marketed device.

Reproducibility was evaluated using control urines containing drug concentrations above and below the stated cut-off. Negative controls were also tested. The results confirmed the reproducibility of the STARTOX™ Drugs of Abuse Screening Test.

510k Submission  
Starplex, Inc.  
Request for additional data  
K022388

**CONCLUSION:**

STARTOX™ Drug of Abuse Screening Test is substantially equivalent to ABMC Rapid Drug Screen multi-drug products (K012159, K002447) previously cleared for market as demonstrated by the results obtained in the studies.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

FEB 27 2003

Starplex Scientific, Inc.  
c/o Ms. Fran White  
Regulatory Consultant  
MDC Associates  
163 Cabot Street  
Beverly, MA 01915

Re: k022388  
Trade/Device Name: Startox Drug of Abuse Screening Test  
Regulation Number: 21 CFR 862.3610  
Regulation Name: Methamphetamine test system  
Regulatory Class: Class II  
Product Code: DJC; DJG; DIO; LDJ  
Dated: December 4, 2002  
Received: December 6, 2002

Dear Ms. White:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of In Vitro Diagnostic Device Evaluation and Safety at (301) 594-3084. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive, slightly slanted style.

Steven I. Gutman, M.D., M.B.A.  
Director  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and  
Radiological Health

Enclosure

510k Submission  
Starplex, Inc.  
Request for additional data  
K022388

**510(K) NUMBER: K022388**

**Device Name: STARTOX Drug of Abuse Screening Test**


**Indication for Use:**

“STARTOX” Drugs of Abuse Screening Test is a one-step lateral flow immunoassay intended for the simultaneous detection of multiple drug analytes in urine. “STARTOX” is intended for use in the qualitative detection of drugs of abuse at the following Substance Abuse Mental Health Services Administration (SAMHSA) recommended levels:

Compound	Abbreviation	Level
Methamphetamine ((+/-)methamphetamine HCl)	METH	1000 ng/ml
Opiates (morphine)	OPIATES	300 ng/ml
Cocaine (benzoylecgonine)	COCAINE	300 ng/ml
Cannabinoids (11-nor- $\Delta^9$ -THC-9-carboxylic-acid)	THC	50 ng/ml

**STARTOX™ Drugs of Abuse Screening Test provide only a preliminary qualitative test result. Use a more specific alternate quantitative analytical method to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Apply clinical and professional judgment to any drug of abuse test result, particularly when preliminary positive results are obtained.**

**For professional use only.**

  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K022388

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

.....

**Concurrence of CDRH, Office of Device Evaluation (ODE)**

Prescription Use ✓  
(Per 21 CFR 801.109)

**OR**

Over The Counter Use \_\_\_\_\_  
(Optional Format 1-2-96)