



MAY 1 5 2003

510(k) Summary:

MorpheusTM 1, Automated Sleep Study Scoring and Data Management System

510(k) Number:

Company Name:

WideMed Ltd.

Contact Person: David Solomon, Ph.D.,

R&D Director

Telephone: +972-8-690-9488.

Fax:

+972-8-690-9489

Trade Name:

MorpheusTM 1, Automated Sleep Study Scoring and Data Management System.

Classification name: Breathing Frequency Monitor

Classification: MNR

Predicate Device:

Compumedics Sleep Monitoring System, Compumedics Sleep Pty. Ltd, Australia cleared under 510(k) no. K955841.

Indications for Use:

The MorpheusTM 1 Automated Sleep Study Scoring and Data Management System is a computer program (software) intended for use as an aid for the diagnosis of sleep and respiratory disorders.

The MorpheusTM 1, Automated Sleep Study Scoring and Data Management System is intended to be used for analysis (automatic scoring and manual rescoring), display, redisplay (retrieve), summarize, reports generation and networking of data received from monitoring devices typically used to evaluate sleep and respiratory related sleep disorders.

This device is to be used under the supervision of a physician.

Substantial Equivalence:

The MorpheusTM 1, Automated Sleep Study Scoring and Data Management System has the same intended use and the same principle of operation as the Compumedics Sleep Monitoring System, cleared under 510(k) no. K955841 and is therefore substantially equivalent to that device.



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

MAY 1 5 2003

Mr. David Solomon WideMed Limited Omer Industrial Park, Building 8c POB 3002 Omer 84965 ISRAEL

Re: K022506

Trade/Device Name: MorpheusTM 1 Automated Sleep Study Scoring and

Data Management System

Regulation Number: 868.2375

Regulation Name: Ventilatory Effort Recorder

Regulatory Class: II Product Code: MNR Dated: March 9, 2003 Received: March 14, 2003

Dear Mr. Solomon:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal</u> Register.

Page 2 – Mr. David Solomon

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4646. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Susan Runner, DDS, MA

Interim Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices

Putine Cuantiffee

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if known):	K022506
Device Name:	Morpheus TM 1, Automated Sleep Study Scoring and Data Management
	System
Indications for Use:	The Morpheus TM 1 Automated Sleep Study Scoring and Data Management System is a computer program (software) intended for use as an aid for the diagnosis of sleep and respiratory disorders.
	The Morpheus TM 1, Automated Sleep Study Scoring and Data Management System is intended to be used for analysis (automatic scoring and manual re-scoring), display, redisplay (retrieve), summarize, reports generation and networking of data received from monitoring devices typically used to evaluate sleep and respiratory related sleep disorders.
	This device is to be used under the supervision of a physician.
IF NEEDED)	WRITE BELOW THIS LINE -CONTINUE ON ANOTHER PAGE
Concurrence of CDRI	H, Office of Device Evaluation (ODE):
	(Division Sign-Off) Division of Anesthesiology, General Hospital, Infection Control, Dental Devices 510(k) Number: 1022506
Prescription Use	OR Over the Counter Use