K022569/A

AUG 2 9 2002

Attachment I 510(K) Summary

Prolite / Plasmalite MPX Pulsed Light System

This 510(K) Summary of safety and effectiveness for the Prolite / Plasmalite MPX Pulsed Light System is submitted in accordance with the requirements of the SMDA 1990 and following guidance concerning the organization and content of a 510(K) summary.

Applicant:

Medical Bio Care AB

Address:

Lona Knapes gata 5 421 32 Vastra Frolunda,

Sweden

Contact Person:

Morgan Gustafsson

Telephone / Fax / Email

46.31.709.30.70 - Phone 46.31.709.30.79 - Fax

morgan@medicalbiocare.com

Preparation Date:

August 11, 2002

Device Trade Name:

Prolite / Plasmalite MPX Pulsed Light Sys

Common Name:

Pulsed Light for Photothermolysis

Classification Name:

Instrument, Surgical, Powered, laser

79-GEX, 21 CFR 878-48

Legally Marketed Predicate Device:

ProLite Pulsed Light System

K number 010928

Description of the Prolite / Plasmalite

MPX Pulsed Light System

The Prolite / Plasmalite MPX Pulsed Light System delivers pulsed light at a wavelength of 550 nanometers. The device consists of three interconnected sections: The cabinet which houses the power supply, the cooling system and the microcontroller, the umbilical to the handpiece, and the

handpiece, which houses the waveguide

Intended use of the Prolite / Plasmalite

MPX Pulsed Light System

The Prolite / Plasmalite MPX Pulsed Light System is

indicated for the treatment of vascular lesions.

Performance Data:

None

Conclusion:

The ProLite / Plasmalite MPX Pulsed Light System is substantially equivalent to other existing pulsed light systems in commercial distribution for treatment of vascular

lesions in Dermatology and Plastic Surgery.



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

AUG 2 9 2002

Medical Bio Care Sweden AB c/o Connie White Hoy 908 Stetson Street Woodland, California 95776

Re: K022569

Trade/Device Name: ProLite/Plasmalite MPX Pulsed Light System

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser Surgical Instrument for Use in General and Plastic

Surgery and in Dermatology

Regulatory Class: Class II

Product Code: GEX Dated: July 25, 2002

Received: August 2, 2002

Dear Ms. Hoy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for <u>in vitro</u> diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sinderely yours

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

INDICATION FOR USE STATEMENT

510(k) Number: K0225	569	
Device Name: ProLite / Plas	smalite MPX Pulsed Light System	
Indications for Use:	• • • • • • • • • • • • • • • • • • •	
The ProLite Pul- treatment of vasc	sed Light System is intended to be used in cular lesions.	the
(Please do not write b	below this line - Continue on another page if needed)	
Concurrence	of CDRH, Office of Device Evaluation (ODE)	
Prescription Use	OR Over-the-Counter Use	
(per 21 CFR 801.109)	(Division Sign-Off) Division of General, Restorative	
	and Neurological Devices	
	510(k) Number <u>KOZZ</u> 56 9	