



July 31<sup>st</sup>, 2001

## **Special 510(k) Summary**

**Echo-View 5.x**

**Easy-View 2.x**

**Omni-View 2.x**

**Cardio-View 1.x**

**LV Analysis 1.x**

**Surgical View 1.x**

OCT 17 2002

### **Name and Address**

TomTec Imaging Systems GmbH  
Edisonstrasse 6  
D-85716 Unterschleissheim

### **Contact Person**

Florian Eisenberger  
Director, Regulatory Affairs & Quality Assurance  
Phone ++49-89-32175-830  
fax ++49-89-32175-750

### **Common, Classification & Proprietary Names**

Common Name:	Digital Ultrasound Image Analysis System
Classification Name:	Ultrasonic Pulsed Echo Imaging System
Proprietary Name(s):	Echo-View 5.x Easy-View 2.x Omni-View 2.x Cardio-View 1.x LV Analysis 1.x Surgical View 1.x

### **Predicate Device**

TomTec Echo-View K993398

### **Device Description**

The Review Software products

- Echo-View 5.x
- Easy-View 2.x
- Omni-View 2.x
- Cardio-View 1.x
- LV Analysis 1.x
- Surgical View 1.x

are software modules for high performance computer systems based on Microsoft Windows 2000/XP™ operating system standards. These Review Software products are proprietary software for the analysis, storage, retrieval and reconstruction of digitized ultrasound B-mode images and Color Doppler images. The data can be acquired by a TomTec acquisition station or by a 3D capable ultrasound system. The result of acquired images allows a 3-dimensional volume to be reconstructed by Echo-View. The digital 3D / 4D data set can be used for 2D and 3D measurements.

### **Intended Use**

The Review Software products

- Echo-View 5.x
- Easy-View 2.x
- Omni-View 2.x
- Cardio-View 1.x
- LV Analysis 1.x
- Surgical View 1.x

are intended to retrieve, analyze and store digital ultrasound images and Color Doppler images for computerized 3-dimensional and 4-dimensional (dynamic 3D) image processing.

Echo-View 5.x , Easy-View 2.x, Omni-View 2.x, Cardio-View 1.x, LV Analysis 1.x and Surgical View 1.x can import certain digital 2D or 3D image file formats for 3D tomographic reconstructions and surface rendering. It is intended as a general purpose digital 3D ultrasound image processing tool for cardiology, radiology, neurology, gastroenterology, urology, surgery, obstetrics and gynecology.

### **Technological Characteristics Comparison**

Echo-View 5.x , Easy-View 2.x, Omni-View 2.x, Cardio-View 1.x, LV Analysis 1.x and Surgical View 1.x are modified version and follow-up products of the filed Echo-View / Easy-View system, which has been transferred to Windows 2000/XP™ operating system standards. The graphic user interface has been improved for faster and easier application.

While Echo-View includes the full functionality the products Cardio-View, LV Analysis, Surgical View, Easy-View and Omni-View are subsets of the Echo-

View Software where the functionality is limited for easier software handling with respect to the software application.

**Test Discussion**

Testing was performed according to internal company procedures. Software testing and validation were done at the module and system level according to written test protocols established before testing was conducted. Test results were reviewed by designated technical professionals before software proceeded to release.

**Test Conclusions**

Test results support the conclusion that actual device performance satisfies the design intent. Actual device performance as tested internally conforms to the system performance specifications.

July 31<sup>st</sup>, 2002



Florian Eisenberger



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

OCT 17 2002

Mr. Florian Eisenberger  
Manager Regulatory Affairs  
TomTec Imaging System  
Edisonstrasse 6  
D-85716 Unterschleissheim  
Germany

Re: K022824

Trade/Device Name: TomTec Echo-View 5.x,  
Easy-View 2.x, Omni-View 2.x, Cardio-View 1.x,  
LV Analysis 1.x, and Surgical View 1.x  
Regulation Number: 21 CFR 892.2050  
Regulation Name: Picture archiving and  
communications system  
Regulation Number: 21 CFR 892.1560  
Regulation Name: Ultrasonic pulsed echo imaging system  
Regulatory Class: II  
Product Code: 90 LLZ and IYO  
Dated: September 16, 2002  
Received: September 24, 2002

Dear Mr. Eisenberger:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

510(k) Number (if known): \_\_\_\_\_

Device Name: TomTec Echo-View 5.x, Easy-View 2.x, Omni-View 2.x  
Cardio-View 1.x, LV Analysis 1.x, Surgical-View 1.x

Indications For Use

The products:


- Echo-View 5.x
- Easy-View 2.x
- Omni View 2.x
- Cardio-View 1.x
- LV Analysis 1.x
- Surgical View 1.x

are intended to retrieve, analyze and store digital ultrasound images and Color Doppler images for computerized 3-dimensional and 4-dimensional (dynamic 3D) image processing.

Echo-View 5.x , Easy-View 2.x, Omni-View 2.x, Cardio-View 1.x, LV Analysis 1.x and Surgical View 1.x can import certain digital 2D or 3D image file formats for 3D tomographic reconstructions and surface rendering. It is intended as a general purpose digital 3D ultrasound image processing tool for cardiology, radiology, neurology, gastro-enterology, urology, surgery, obstetrics and gynecology.

(PLEASE DO NOT WRITE BELOW LINE LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number 1K022824

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_