

NOV 4 2002

1022948

Section Seven: 510(k) Summary for Vaginal Speculums

Submitted by: Medical Action Industries, Inc. (MAI)
25 Heywood Road
Arden, North Carolina 28704
Telephone: 828-681-8820

Contact Person: Mr. Dennis Kanupp

Date Prepared: September 4, 2002

Proprietary Name: Not Applicable

Common Name: Disposable vaginal speculum

Classification Name: Speculum, Vaginal, Non-metal

Predicate Device: Leespec Disposable Vaginal Speculum

Description of Device: The MAI vaginal speculum is a non-metal (polystyrene), hand held device used to expose the interior of the vagina.

Intended use of the Device: To be used by medical professionals to expose the interior of the vagina to facilitate visualization during gynecological and obstetrical procedures.

Technical Characteristics: The MAI vaginal speculum has the same technological characteristics as and is substantially equivalent to the Leespec Disposable Vaginal Speculum (K021017), manufactured by ITL Corporation.



NOV 4 2002

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Medical Action Industries, Inc.
% Ms. Sharon Suess Graham
Consultant
Ryan Dietrich Engineering
2 Summit Drive
ARDEN NC 28704

Re: K022948
Trade/Device Name: Vaginal Speculum
Regulation Number: 21 CFR 884.4530
Regulation Name: Obstetric-gynecologic specialized
manual instrument
Regulatory Class: II
Product Code: 85 HIB
Dated: September 4, 2002
Received: September 5, 2002

Dear Ms. Graham:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Section Six: Indications for Use Statement

510(k) Number: K 022948

Device Name: Vaginal Speculum

Indications for Use: Vaginal Speculums are single use and shall be used to expose the interior of the vagina.


PLEASE DO NOT WRITE BELOW THIS LINE- CONTINUE ON ANOTHER PAGE IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The Counter Use ☐


(Division Sign-Off)
Division of Reproductive, Abdominal,
and Radiological Devices
510(k) Number K022948