

NOV 22 2002

K023636

Section 807.87 (h) A 510(k) Summary as described in Section 807.92 or a 510(k) statement as described in 807.93

Premarket Notification [510(k)] Summary as required by 21 CFR 807.92

Date summary was prepared:

October 29, 2002

Submitter's Name:

Southeastern Radiation Products, Inc.
2651 N. Design Ct.
Sanford, FL 32773

Contact Person:

Richard L. Sweat
President
Phone: 407-330-3300
FAX: 407-322-7546
Email: rsweat@seradiation.com

Device Name:

Blocking Tray

Classification Name:

Accessory to a Linac, Class II

Predicate Device(s):

Accessory to Various Linac Manufacturers including:

Varian	K913174
	K913119

K023636

Siemens	K982513
Philips (Elekta)	K963624
GE/CGR	K833941
Mitsubishi	K970137
Mick-Radio	K991107
Med-Tec Inc	K944315

Intended Use:

Blocking Trays are to be used as an accessory to a Linac. Custom blocks are mounted onto Blocking trays and slid into the Linac into their accessory holder thereby stenciling the radiation beam onto the patient. Compensators are also mounted on blocking trays and slid into the Linac changing the radiation beam fluence unique to a patient's anatomy.

Technological Characteristics:

The technological characteristics of blocking trays are discussed in Southeastern Radiation Products's attached documents:

DMR-02	Device Master Record – Custom Blocking Trays
PD-02	Processing a Blocking Tray Work Order
Acrylic MSDS	
Polycarbonate MSDS	



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

NOV 22 2002

Mr. Richard L. Sweat
President
Southeastern Radiation Products, Inc.
2651 N. Design Court
SANFORD FL 32773

Re: K023636
Trade/Device Name: Blocking Trays
Regulation Number: 21 CFR 892.5050
Regulation Name: Medical charged-particle
radiation therapy system
Regulatory Class: II
Product Code: 90 IYE
Dated: October 29, 2002
Received: October 30, 2002

Dear Mr. Sweat:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon
Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

510 (k) Number (if known): K023636

Device Name: Blocking Trays

Indications for Use:

This device is intended for mounting blocks and/or compensators/modulators for insertion into a Linear Accelerator head.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

Nancy C. Burdick
(Division Sign-Off)
Division of Reproductive, Abdominal,
and Radiological Devices
510(k) Number K023636