



MAY 30 2003

K023887

Page 1 of 2

## **510(k) SUMMARY**

### **Inion OTPS™ Biodegradable Mini Plating System**

#### **MANUFACTURER**

Inion Ltd.  
Lääkärintäti 2  
FIN-33520 Tampere

Contact Person:  
Hanna Marttila  
Regulatory Affairs Manager  
Phone: +358 3 2306 600  
Fax: +358 3 2306 601  
[Hanna.Marttila@Inion.fi](mailto:Hanna.Marttila@Inion.fi)

#### **DEVICE NAME**

Trade name: Inion OTPS™ Biodegradable Mini Plating System  
Common/Usual Name: Bone Plating System  
Classification Name: Bone Plate, Screw

#### **ESTABLISHMENT REGISTRATION NUMBER**

9710629

#### **DEVICE CLASSIFICATION AND PRODUCT CODE**

Classification panel: Orthopedic  
Regulatory Class: Class II  
21 CFR 888.3040 – Smooth or threaded metallic bone fixation fastener, 87 HWC (screw)  
21 CFR 888.3030 – Single multiple component metallic bone fixation appliances and accessories, 87 HRS (plate)

#### **PREDICATE DEVICES**

- (1) Biomet Inc.; ReUnite Bone Screw (K992301)
- (2) Biomet Inc.; LactoSorb Hand System (K991763)
- (3) Syntehes(USA); Synthes Compact Hand 1.0/1.3 (pre-amendment device)

## **DEVICE DESCRIPTION AND PRINCIPLES OF OPERATION**

The Inion OTPS™ Biodegradable Mini Plating System is intended for use in fixation of non-comminuted diaphyseal fractures of the metacarpal, proximal phalangeal, middle phalangeal and osteotomies in the presence of appropriate immobilization.

The Inion OTPS™ Biodegradable Mini Plating System is consisted of biodegradable screws and plates intended to be used in fractures of phalanges and metacarpal. Implants will be utilized by using common surgical techniques. The Inion OTPS™ Biodegradable Mini Plating System is made of resorbable polylactic acid / trimethylenecarbonate copolymers [Poly (L-lactide-co-D,L-lactide) and poly (L-lactide-co-trimethylenecarbonate)]. The Inion OTPS™ Biodegradable Mini Plating System implants gradually lose their strength during 18-36 weeks in vivo with complete strength loss and resorption within two to four years.

## **EQUIVALENCE TO MARKETED PRODUCTS**

The Inion OTPS™ Biodegradable Mini Plating System is substantially equivalent to biodegradable implants, intended for metacarpal and phalangeal fracture fixation procedures, which have received 510(k) clearance. Inion OTPS™ Biodegradable Mini Plating System LactoSorb Hand System (K991763) and ReUnite Screw (K992301) have the same intended use and principles of operation and very similar design characteristics. Mechanical testing demonstrates that the device is substantially equivalent to the predicate ones. The differences between the Inion OTPS™ Biodegradable Mini Plating System and predicate devices do not raise new questions of safety and effectiveness



MAY 3 0 2003

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Hanna Marttila  
Regulatory Affairs Manager  
Inion Ltd.  
Lääkärinkatu 2  
Fin-33520 Tampere  
Finland

Re: K023887

Trade/Device Name: Inion OTPS™ Biodegradable Mini Plating System

Regulation Number: 21 CFR 888.3040 and 888.3030

Regulation Name: Smooth or threaded metallic bone fixation fastener, and Single/multiple component metallic bone fixation appliances and accessories

Regulatory Class: II

Product Code: HWC and HRS

Dated: March 4, 2003

Received: March 6, 2003

Dear Ms. Marttila:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97) you may obtain.

Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address

<http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "for Celia M. Witten". The signature is written in a cursive, flowing style.

Celia M. Witten, Ph.D., M.D.  
Director

Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**D STATEMENT OF INDICATIONS FOR USE****Applicant:** Inion Ltd.**510(k) Number:** K023887**Device Name:** Inion OTPS™ Biodegradable Mini Plating System**Indications:**

The Inion OTPS™ Biodegradable Mini Plating System is intended for use in fixation of non-comminuted diaphyseal fractures of the metacarpal, proximal phalangeal, middle phalangeal and osteotomies in the presence of appropriate immobilization.

**Contraindications:**

Inion OTPS™ Biodegradable Mini Plating System is not intended for use in and is contraindicated for:

1. Active or potential infection
2. Patient conditions including limited blood supply, insufficient quantity or quality of bone; and where patient cooperation cannot be guaranteed. (e.g., alcoholism, drug abuse)
3. Load bearing procedures.

Prescription use yes Over the Counter use No

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Per 21 CFR 801.109)

for Mark N. Miller  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

510(k) Number K023887

Date: 15.11.2002  
Status: Final