

510(k) Summary – K024348

SUBMITTER: DePuy AcroMed, Inc.
325 Paramount Drive
Raynham, MA 02767

CONTACT PERSON: Lisa A. Gilman

DATE PREPARED: December 16, 2002

CLASSIFICATION NAME: Appliance, Fixation, Spinal Interlaminar
Orthosis, Spinal Pedicle Fixation

PROPRIETARY NAME: MONARCH Spine System

PREDICATE DEVICES: MONARCH Spine System (K021335, K010576)

DEVICE DESCRIPTION: Addition of 5.5mm diameter rod, screws, hooks and
various connectors to the MONARCH Spine System.

The MONARCH Spine System also contains Class 1
manual surgical instruments and cases that are
considered exempt from premarket notification.

INTENDED USE: The MONARCH Spine System is a pedicle screw
system intended to provide immobilization and
stabilization of spinal segments in skeletally mature
patients as an adjunct to fusion in the treatment of the
following acute and chronic instabilities or deformities
of the thoracic, lumbar, and sacral spine:
degenerative spondylolisthesis with objective
evidence of neurological impairment, fracture,
dislocation, scoliosis, kyphosis, spinal tumor, and
failed previous fusion (pseudarthrosis).

The MONARCH Spine System is also indicated for
pedicle screw fixation for the treatment of severe
spondylolisthesis (Grades 3 and 4) of the L5-S1
vertebra in skeletally mature patients receiving fusion
by autogenous bone graft having implants attached to
the lumbar and sacral spine (L3 to sacrum) with
removal of the implants after the attainment of a solid
fusion.

The MONARCH Spine System is also a hook and sacral/ilic screw fixation system of the noncervical spine indicated for degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies), spondylolisthesis, trauma (fracture and/or dislocation), spinal stenosis, deformities (scoliosis, lordosis and/or kyphosis), tumor, and previous failed fusion (pseudarthrosis).

The MONARCH Spine System when used with pedicle screws is indicated for degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies). Levels of fixation are for the thoracic, lumbar and sacral spine.

MATERIALS:

Manufactured from ASTM F-136 implant grade titanium alloy.

**PERFORMANCE
DATA:**

Performance data were submitted to characterize the MONARCH Spine System.



JAN 28 2003

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Frank Maas
Director, Regulatory Affairs
Depuy AcroMed
325 Paramount Drive
Raynham, MA 02767-0350

Re: K024348
Trade Name: MONARCH Spine System – Addition of 5.5 mm Components
Regulation Number: 21 CFR 888.3070 (b)(1) and 888.3050
Regulation Name: Pedicle screw spinal system, and spinal Interlaminar fixation orthosis
Regulatory Class: III
Product Code: MNH, MNI, KWP
Dated: December 26, 2002
Received: December 30, 2002

Dear Mr. Maas:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", with a stylized flourish at the end.

Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K024348

Device Name: MONARCH Spine System

Indications For Use:

The MONARCH Spine System is a pedicle screw system intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of the thoracic, lumbar, and sacral spine: degenerative spondylolisthesis with objective evidence of neurological impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, and failed previous fusion (pseudarthrosis).

The MONARCH Spine System is also indicated for pedicle screw fixation for the treatment of severe spondylolisthesis (Grades 3 and 4) of the L5-S1 vertebra in skeletally mature patients receiving fusion by autogenous bone graft having implants attached to the lumbar and sacral spine (L3 to sacrum) with removal of the implants after the attainment of a solid fusion.


The MONARCH Spine System is also a hook and sacral/iliac screw fixation system of the noncervical spine indicated for degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies), spondylolisthesis, trauma (fracture and/or dislocation), spinal stenosis, deformities (scoliosis, lordosis and/or kyphosis), tumor, and previous failed fusion (pseudarthrosis).

The MONARCH Spine System when used with pedicle screws is indicated for degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies). Levels of fixation are for the thoracic, lumbar and sacral spine.

(Please do not write below this line - continue on another page if needed)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use: _____ OR Over-The-Counter Use: _____
(Per 21 CFR 801.109)


Division Sign-Off)
Restorative
and Medical Devices

DePuy AcroMed, Inc.
Special 510K

510(k) Number K024348 4