

K030610

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS

I. **Submitter:** Tim Ehrecke, Portable Ophthalmic Devices, Inc.
4413 Winston Place,
Bettendorf, Iowa 52722.
Ph: 563-529-3536 FAX 563-449-9184

MAR 21 2003

II. **Classification Names and Numbers:** Ultrasonic Pulsed Echo Imaging System,
90-IYO, 21CFR892.1560

III. **Common/Usual Name:** Ultrasonic Corneal Pachymeter

IV. **Proprietary Names:** POD Pac-20m™

V. **Establishment Registration Number:** in progress

VI. **Classification:** Class II, Tier II. Described in 21CFR892.1560

VII. **Substantial Equivalence:** We believe the POD Pac-20m™ is substantially equivalent to the classified device described in 21CFR892.1560, "Ultrasonic Pulsed Echo Imaging System," and to other ultrasound systems that have been cleared by the 510(k) process, such as the Quantel "Pocket" Ultrasonic Pachymeter (K993674).

The 510(k) Substantial Equivalence Decision-making Process (detailed) from ODE Guidance Memorandum #86-3 was followed as described below:

1. These products have the same intended use, to measure the thickness of the cornea.
2. The technological characteristics of this device are the same as those for the predicate device except for the models of some of the components, and the user interface.
3. The materials from which the patient-contact portions of the device are made are the same as in the predicate device.
4. The acoustic output of this device is similar to that of predicate devices and below the preamendment levels described in the guidance document "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers."



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAR 21 2003

Mr. Tim Ehrecke
President and CEO
Portable Ophthalmic Devices, Inc.
4413 Winston Place
BETTENDORF IA 52772

Re: K030610

Trade Name: POD Pac-20m Corneal Pachymeter System
Regulation Number: 21 CFR 892.1560
Regulation Name: Ultrasonic pulsed echo imaging system
Regulatory Class: II
Product Code: 90 IYO
Dated: February 25, 2003
Received: February 26, 2003

Dear Mr. Ehrecke:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the POD Pac-20m Corneal Pachymeter System, as described in your premarket notification:

Transducer Model Number

Blatek AT15399/AT15387 (20 MHz)

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:


Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center (HFZ-401)
9200 Corporate Boulevard
Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801, please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

If you have any questions regarding the content of this letter, please contact Rodrigo C. Perez at (301) 594-1212.

Sincerely yours,



Nancy C. Brogdon
Director, Division of Reproductive,
Abdominal and Radiological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure(s)

Diagnostic Ultrasound Indications for Use Form

510(k) Number: K030610

Device Name: POD Pac-20m Corneal Pachymeter System

Intended Use: To measure thickness of the cornea in the human eye using ultrasound energy.

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)	Other (specify)
Ophthalmic	N									
Fetal										
Abdominal										
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (specify)										
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Peripheral Vascular										
Laparoscopic										
Musculoskeletal Conventional										
Musculoskeletal Superficial										
Other (specify)										

N= new indication; P= previously cleared by FDA; E= added under Appendix E

Additional Comments:

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109) ✓

(Division Sign-Off)

Division of Reproductive, Abdominal,
and Radiological Devices

510(k) Number

K030610