

MAY 28 2003

K030667 p1/1

**510(k) Summary of Safety and Effectiveness**

- (1) **Submitter's name:** Encore Medical, L.P.  
**Submitter's address:** 9800 Metric Blvd, Austin, TX 78758  
**Submitter's telephone number:** (512) 834-6255  
**Contact person:** Debbie De Los Santos  
**Date summary prepared:** February 19, 2003
- (2) **Trade or proprietary device name:** Encore Anterior Spinal Staple System  
**Common or usual name:** Anterior Spinal Staple System  
**Classification name:** 888.3060 – Spinal Intervertebral Body Fixation Orthosis
- (3) **Predicate devices:** AcroMed KASS (K974757)  
Blackstone Medical (K022605)
- (4) **Subject device description:**  
The anterior vertebral body staple is designed to work in conjunction with the any existing spinal fixation systems utilizing 6.0mm to 6.5mm titanium screws. The staples are manufactured from Titanium alloy and are available as a single or double staple. The single staple has one hole for screw fixation with a 6.0mm – 6.5mm titanium screws. The double staples are designed in pairs with specific superior and inferior components. Each has two screw holes allowing screws to pass through at a converging angle.
- (5) **Subject device intended use:**  
The Encore Anterior Spinal Staple System is intended for anterior vertebral body fixation to the T3-L5 spine. This system is indicated to be used with an anterior spinal system that consists of screws that are between 6.0 and 6.5mm diameter manufactured out of Ti alloy.
- (6) **Basis for Substantial Equivalence:**  
The Encore Anterior Spinal Staple System is similar in design, indications and materials to AcroMed's Kaneda Anterior Spinal System (K974757) and Blackstone Spinal Fixation System (K022605).

MAY 28 2003  
DH



MAY 28 2003

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Debbie De Los Santos  
Supervisor, Regulatory/QA Manager  
Encore Orthopedics  
9800 Metric Boulevard  
Austin, Texas 78758

Re: K030667  
Trade Name: Encore Anterior Spinal Staple System  
Regulatory Number: 21 CFR 888.3060  
Regulation Name: Spinal intervertebral body fixation orthosis  
Regulatory Class: II  
Product Code: KWQ  
Dated: February 28, 2003  
Received: March 3, 2003

Dear Ms. De Los Santos:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Ms. Debbie De Los Santos

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "for Celia M. Witten". The signature is written in a cursive, flowing style.

Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known): K030667

Device Name: Encore Anterior Spinal Staple System

Indications For Use:

**Encore Anterior Spinal Staple System**  
**Indications For Use**

The Encore Anterior Spinal Staple System is intended for anterior vertebral body fixation to the T3-L5 spine. This system is indicated to be used with an anterior spinal system that consists of screws that are between 6.0 and 6.5mm diameter manufactured out of Ti alloy.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)


Concurrence of CDRH, Office of Device Evaluation (ODE)

---

Prescription Use \_\_\_\_\_  
(per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_  
(Optional Format 1-2-96)\_

  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

510(k) Number K030667