#### Appendix I, Page 1 (added)

### 510(k) SUMMARY FOR

#### K030993

# NON-STERILE GREEN COLORED POWDERED LATEX PATIENT EXAMINATION GLOVES WITH & WITHOUT SOUR APPLE/PEPPERMINT SCENTS AND WITH A PROTEIN LABELING CLAIM (<200ug/g)

Confains 200 migm onless of total Water Extractable Protein pergram.

Submitted for:

SGMP COMPANY, LTD.

Submitted by:

**TUCKER & ASSOCIATES** 

Official Correspondent for SGMP COMPANY, LTD.

JANNA P. TUCKER, President - CEO

198 Avenue de la D'emerald Sparks, NV 89434-9550

Phone:

775-342-2612

Fax:

775-342-2613

E-Mail:

Tuckerjan@aol.com

This device is substantially equivalent to K000671, which is another of SGMP's Colored, and/or scented powdered latex gloves with protein labeling (<200ug/g).

Revised 06-04-03

# 2. Physical Properties (ASTM-D3578-01aE2 Standard Specification for Latex Exam Gloves)

LOT#		LE STREN	IGTH	ULTIMATE ELONGATION				
0238	UNAGED		AGED		UNAGED		AGED	
TESTED X-SMALL	<u>SGMP</u>	<u>ASTM</u>	SGMP	<u>ASTM</u>	SGMP	<u>ASTM</u>	<u>SGMP</u>	<u>ASTM</u>
	22.6	14.0	20.8	14.0	820	700	950	500
SMALL	21.9	14.0	22.5	14.0	920	700	920	500
MEDIUM	21.7	14.0	22.8	14.0	970	700	880	500
LARGE	21.3	14.0	20.7	14.0	950	700	896	500

#### 3. Water Tight Test Data

BATCH NUMBER	DATE TESTED	SAMPLING SIZE	LEAK STATUS	NUMBER LEAKED
Unaged Smpl				
0238 XS	10 Feb 03	125	No	0
S		125	Yes	1
M		125	Yes	2
L		125	Yes	1
Aged Smpl				
0238 XS	18 Feb 03	125	No	0
S		125	Yes	1
M		125	Yes	2
L		125	Yes	1

The above figures are within the ASTM D-3578-01aE2 requirements for latex exam gloves of 2.5% AQL.

#### 4. Biocompatibility

**BIOCOMPATIBILITY TESTS** 

Test results indicate that the gloves passed the biocompatibility tests for gloves.

#### 5. Residual Protein Level

TESTS	FDA ALLOWABLE LEVEL	CLAIMED LEVEL
ASTM D 5712-99	-	< 200 μg/ <b>3</b> Range: 91 –123 <u>μg/9</u>
		Mean: 112 μg/g

Revised 5.28.03



JUN 1 2 2003

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

SGMP Company, Limited C/O Ms. Janna P. Tucker Tucker & Associates 198 Avenue de la D' emerald Sparks, Nevada 89434-9550

Re: K030993

Trade/Device Name: Non-Sterile Powdered Green Colored Powdered Latex Patient with or without Sour Apple, Peppermint Scents, Contains 200 Micrograms or Less

of Total Extractable Protein Per Gram

Regulation Number: 880.6250

Regulation Name: Patient Examination Gloves

Regulatory Class: I Product Code: LYY Dated: May 28, 2003 Received: May 29, 2003

#### Dear Ms. Tucker:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal</u> Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4618. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Susan Runner, DDS, MA

Interim Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

## INDICATIONS FOR USE

A	D	D	T	T		A	N	T	•
$\boldsymbol{A}$		r	L	11	v.	$\sim$	JΥ	ı	•

SGMP COMPANY, LTD.

510(k) NUMBER:

**DEVICE NAME:** 

NON-STERILE GREEN COLORED

POWDERED LATEX PATIENT

**EXAMINATION GLOVES WITH &** 

WITHOUT SOUR

APPLE/PEPPERMINT SCENTS AND

WITH A PROTEIN LABELING CLAIM

(<200 uG/G) Contains 200 migm, of to fail water extractable protein per gram

A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDEDI

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Anesthesiology, General Hospital, Infection Control, Dental Devices

510(k) Number:\_

Over-The-Counter Use

Prescription Use\_ (Per 21 CFR 801,109)

(Optional Format 1-2-96)

Kevised 06-04-03