

K030996

## Section E - 510(k) Summary

**Submitted by:**

MedSphere International, Inc.

48511 Warm Springs Blvd., Suite 212, Fremont, CA 94539

Tel: 510-656-8232 Fax: 510-656-8816

**Contact Person:**

Eric Kao, Vice President of Quality and Regulatory Affairs

**Date Summary Prepared:**

March 28, 2003

**Name of the Device:**

Trade / Proprietary Name: MSI SA Electrodes

Common / usual Name: Electrosurgical electrode

Classification Name: Electrosurgical cutting and coagulation device (21CFR878.4400)

**Predicate Devices:**

MSI TA Electrodes

**Description:**

The MSI SA series of disposable electrodes are electrosurgical devices used for coagulation of soft tissue. These devices are designed for percutaneous, laparoscopic or intraoperative use. The device consists of a needle electrode and a handle. For temperature measurement model, the needle contains a thermocouple for monitoring tissue temperature. The needle electrode is housed in an insulated cannula. The handle incorporates an electrical connector for interconnecting to an electrosurgical generator.

**Statement of Intended Use:**

MSI SA Electrode is intended for coagulation of soft tissue

This device is intended for use by qualified medical personnel trained in the use of electrosurgery.

**Comparison to Predicate Devices:**

MSI SA Electrode has been compared to previously 510(k) cleared device with respect to intended use and technological characteristics. Performance testing was done to validate its intended use. The comparison and performance testing results in this 510(k) notification shows MSI SA Electrode is substantially equivalent to predicate device and is safe and effective in its intended use.



MAY 30 2003

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Eric Kao  
Vice President of Quality  
and Regulatory Affairs  
MedSphere International, Inc.  
48511 Warm Springs Boulevard, Suite 212  
Fremont, California 94539

Re: K030996  
Trade/Device Name: MSI SA Electrodes  
Regulation Number: 21 CFR 878.4400  
Regulation Name: Electrosurgical cutting and coagulation device and accessories  
Regulatory Class: II  
Product Code: GEI  
Dated: March 28, 2003  
Received: March 31, 2003

Dear Mr. Kao:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Mr. Eric Kao

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97) you may obtain. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,


for Mark N. Milken

Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Section D - Statement of Indications for Use

K030996Indications for Use

Applicant:	MedSphere International, Inc. 48511 Warm Springs Blvd., Suite 212, Fremont, CA 94539
510(k) Number (if known):	
Device Name:	MSI SA Electrodes
Indications For Use:	Indicated for coagulation of soft tissue  These devices are intended for use by qualified medical personnel trained in the use of electrosurgery.
(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE, IF NEEDED)	
(Concurrence of CDRH, Office of Device Evaluation (ODE))	
 (Division Sign-Off) Division of General, Restorative and Neurological Devices 510(k) Number <u>K030996</u>	
Prescription Use <input checked="" type="checkbox"/>	OR Over-the-Counter Use _____ (Per 21 CFR 801.109)
(Optional format 1-2-06)	