

510(K) Summary

Submitter: Cynosure, Inc.
10 Elizabeth Drive
Chelmsford, MA 01824
JUL 08 2003

Contact: George Cho
Senior Vice President of Medical Technology

Date Summary Prepared: April 18, 2003

Device Trade Name: Photolight PL

Common Name: Medical Laser System

Classification Name: Instrument, surgical, powered, laser
79-GEX
21 CFR 878.48

Equivalent Device: EsteLux pulsed light system

Device Description: Photolight PL is a intense pulsed light system, having a Xenon flashlamp located in the handpieces. It is a light source with a range approximately 400 - 1400 nm wavelength.

Emission activation is by footswitch. Overall weight of the laser is 18 Kg, and the size is 20x48x54 cm (HxWxD).

Electrical requirement is 110 VAC, 15A, 50-60 Hz, single phase.

Intended Use: The Photolight PL is indicated for permanent hair reduction, photocoagulation of vascular lesions, photothermolysis of blood vessels, treatment of benign pigmented lesions.

The CoolHand cooling device is intended to provide pre-cooling of the epidermis, to reduce thermal injury to the epidermis, and to reduce patient pain and discomfort associated with light applications.

Comparison: The Photolight PL system has an equivalent indication for uses, the same principle of operation, and essentially the same wavelength range and pulse energy range as the predicate device.

Nonclinical Performance Data: none

Clinical Performance Data: none

Conclusion: The Photolight PL is another safe and effective device for permanent hair reduction, photocoagulation of vascular lesions, photothermolysis of blood vessels, treatment of benign pigmented lesions.

Additional Information: none



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUL 08 2003

Mr. George Cho
Senior Vice President of Medical Technology
Cynosure, Inc
10 Elizabeth Drive
Chelmsford, Massachusetts 01824

Re: K031258
Trade/Device Name: Cynosure Photolight PL
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser surgical instrument for use in general
and plastic surgery and in dermatology
Regulatory Class: II
Product Code: GEX
Dated: June 12, 2003
Received: June 13, 2003

Dear Mr. Cho:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97) you may obtain. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,


for Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K031258

Device Name: Cynosure Photolight PL

Indications For Use:

The Photolight Pulsed Light device is indicated for permanent hair reduction. It is also indicated for photocoagulation of dermatological vascular lesions, photothermolysis of blood vessels (treatment of facial and leg veins), and the treatment of benign pigmented lesions.

The CoolHand cooling device is intended to provide pre-cooling of the epidermis, to reduce thermal injury to the epidermis, and to reduce patient pain and discomfort associated with light applications.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Miriam C. Provost
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices

510(k) Number K031258

Prescription Use ✓

OR

Over-The-Counter Use _____