

XIII. SUMMARY OF SAFETY AND EFFECTIVENESS

510(k) SUMMARY OF SAFETY AND EFFECTIVENESS ESTEEM STERILE POLYISOPRENE SURGICAL GLOVES WITH COATING

Regulatory Affairs Contact:

Erica Sethi

Cardinal Health

1500 Waukegan Road, MP-WM

McGaw Park, IL 60085

Telephone:

(847) 785-3337

Date Summary Prepared:

4/7/03

Product Trade Name:

Undetermined

Common Name:

Surgical Glove

Classification:

Glove, Surgeon's

Predicate Devices:

Esteem Sterile Polyisoprene Surgical Gloves

Description:

Esteem Sterile Polyisoprene Surgical Gloves with Coating are formulated using Synthetic Rubber Latex. These are offered

powder-free and sterile.

Intended Use:

Esteem Sterile Polyisoprene Surgical Gloves with Coating are intended for use in environments within hospitals and other healthcare facilities. The gloves are appropriate for use during invasive and non-invasive medical procedures are similar.

invasive and non-invasive medical procedures requiring sterility. They are intended to be worn by operating room personnel to protect a surgical wound from contamination.

Page 2

Substantial Equivalence:

Esteem Sterile Polyisoprene Surgical Gloves with Coating are substantially equivalent to Esteem Sterile Polyisoprene Surgical Gloves in that they provide the following characteristics:

- same intended use

- same sizes, product features, packaging

- both made of Synthetic Rubber Latex using similar manufacturing

process

Summary of Testing:

<u>Test</u> <u>Result</u>

Primary Skin Irritation Gloves are non-irritating.

Guinea Pig Maximization Gloves do not display any potential for sensitization.

Ultimate Elongation Gloves exceed requirements for rubber surgical gloves & Tensile Strength per ASTM D3577-01.

Barrier Defects Gloves exceed requirements per 21 CFR \$800.20 and ASTM D3577-01, AQL = 1.5.

Data/Test Method Gloves meet powder level requirements for "Powder Free"

designation using ASTM Standard D6124-00-Standard test method for residual powder on medical gloves. Results generated values

below 2 mg of residual powder per glove.



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

OCT 1 5 2003

Ms. Erica Sethi Manager, Regulatory Affairs Cardinal Heath Medical Product and Services 1500 Waukegan Road Building WM Mcgaw Park, Illinois 60085

Re: K031301

Trade/Device Name: Esteem Sterile Polyisoprene Powder Free Surgical Gloves with

Coating

Regulation Number: 878.4460 Regulation Name: Surgeon's Glove

Regulatory Class: I Product Code: KGO Dated: July 23, 2003 Received: July 24, 2003

Dear Ms. Sethi:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4618. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure



Applicant:

1500 Waukegan Road McGaw Park, Illinois 60085-6787 847.473.1500 FAX: 847.785.2460

Cardinal Health

Page 1 of 1

510(k) Number:	K031301
Device Name:	Esteem Sterile Polyisoprene Powder Free Surgical Gloves With Coating
Indications For Use:	These gloves are intended to be worn by operating room personnel to protect a surgical wound from contamination in environments within hospitals and other healthcare facilities. The gloves are appropriate for use during invasive as well as non-invasive medical procedures requiring sterility.
·	•
(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE	
Concurrence of CDRH, Office of Device Evaluation (ODE)	
Ďi	Over-The Counter Use Salue Af Anter Bunch Chief, 10/15/03 ivision Sign-Off) vision of Anesthesiology, General Hospital, ection Control, Dental Devices
51	0(k) Number: <u> </u>