

JUN 17 2003

K031326

Premarket Notification 510(k)

Imagine Reflex

**5. 510 (k) Summary**

Submitter of 510(k):      Wieland Dental + Technik GmbH & Co. KG  
Schwenninger Str. 13  
D-75179 Pforzheim  
Germany  
Phone: +49-7231-3705-0

Contact person:          Dr. Gerhard Polzer  
Phone:                      +49-7231-3705-219  
Fax:                         +49-7231-357959  
e-mail:                     gerhard.polzer@wieland-dental.de

Date of Summary:        2003-04-02

Trade name:              IMAGINE® Reflex

Classification name:     Porcelain powder for clinical use  
Product code:            EIH  
C.D.R section:           872.6660  
Classification:           Class II

Legally marketed  
equivalent device:       Duceram-Metal-Ceramic Dental Porcelain System

510(k) number:          K871808



JUN 17 2003

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Gerhard Polzer, Ph.D.  
Director, Regulatory Affairs  
Wieland Dental + Technik GmbH & Co. KG  
Schwenninger Strasse 13  
D-75179 Pforzheim  
GERMANY

Re: K031326  
Trade/Device Name: Imagine® Reflex®  
Regulation Number: 21 CFR 872.6660  
Regulation Name: Porcelain Powder for Clinical Use  
Regulatory Class: II  
Product Codes: EIH  
Dated: April 08, 2003  
Received: April 25, 2003

Dear Dr. Polzer:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97) you may obtain. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Susan Runner". The signature is fluid and cursive, with a large initial "S" and "R".

Susan Runner, DDS, MA

Interim Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known): K031326

**IMAGINE® Reflex®**

Device Name: \_\_\_\_\_

Indications For Use:

Imagine Reflex is a dental ceramic that can be used by dental technicians to fabricate porcelain-fused-to-metal restorations. It is suitable for veneering dental alloys with a coefficient of thermal expansion [CTE(25 - 500°C):  $13.8 - 15.1 \times 10^{-6} K^{-1}$ ] and a solidus temperature that has to exceed 1000°C (1830°F) as well as electroformed galvanogold.

Imagine Reflex can be used for veneering inlays, partial crowns, single crowns, and bridges.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_

Ran Muly Sa HSR  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital,  
Infection Control, Dental Devices

(Optional Format 1-2-96)

510(k) Number: K 031326