AUG 1 4 2003

ENCLOSURE L

403/8/4 Premarket Notification (510(k) - Mercury Medical StatCO₂meter™ End Tidal CO₂ Detector with Integrated Airway Pressure Manometer.

Non-Confidential Summary of Safety and Effectiveness

Page 1 of 2 June 10, 2003

Mercury Enterprises, Inc./Mercury Medical 11300 49th St. N. Clearwater, Fl. 33762

Tel: (800) 237-6418 Fax: (727) 572-4501

Official Contact:

Wayne Glover

OA/RA Engineer

Proprietary or Trade Name: StatCO₂meter™

Common/Usual Name: End Tidal CO₂ detector with Integrated Airway Pressure

Manometer

Classification: Class II, CCK, 21 CFR 868.1400

Classification Name: Analyzer, Gas, Carbon Dioxide, Gaseous Phase

Device: StatCO₂meter™

Predicate Devices: StatCO₂™ (K021576) and Airway Pressure

Manometer (K954486)

Device Description:

The Mercury Medical StatCO₂meter™ is a fast, durable colorimetric breath indicator for visualization of exhaled CO2 and also measures Airway and PEEP pressure. The StatCO₂meter™ is designed to connect between an endotracheal tube and a breathing device to help verify proper intubation. Exhaled gas passes through the indicator to detect approximate ranges of end-tidal CO₂ by color comparison. The detector may be used during patient transport or in locations where intubations are performed.

ENCLOSURE L

Premarket Notification (510(k) - Mercury Medical StatCO₂meter™ End Tidal CO₂ Detector with Integrated Airway Pressure Manometer.

Non-Confidential Summary of Safety and Effectiveness

Page 2 of 2 June 10, 2003

In	ten	de	ı.	TI	80
ш	LCH	ut	·u	u	SC.

Indicated Use: The Mercury Medical StatCO₂™ End Tidal CO₂ Detector with Integrated Airway Pressure Manometer provides a semi-quantitative visualization of the CO₂ in the patient airway while measuring Airway and PEEP pressure. It is an adjunct in patient assessment, to be used in conjunction with other methods to determine clinical signs and symptoms by or on the order of a physician.

Technical Characteristics: The device has the same technical characteristics as the combination of the following Mercury Medical predicate devices:

StatCO₂™ End Tidal CO₂ Detector Airway Pressure Manometer

Non-Clinical Data: Performance and specifications of the device are consistent with all requirements for this device type specified by ISO 5356-1: 1987 – Anesthetic and Respiratory Equipment-Conical connectors-Part 1: Cones and Sockets. ASTM F1054 – Standard Specification for Conical Fittings of 15mm and 22mm sizes.

Environment of Use: 'Hospital/Transport'

Conclusions: The comparison to the predicate devices demonstrates that the proposed device is safe and effective and is substantially equivalent to the predicate devices.



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

AUG 1 4 2003

Mr. Wayne Glover QA/RA Engineer Mercury Medical 11300-49th Street North Clearwater, Florida 33762-4800

Re: K031814

Trade/Device Name: StatCO₂meterTM

Regulation Number: 868.1400

Regulation Name: Carbon Dioxide Gas Analyzer

Regulatory Class: II Product Code: CCK, CAP Dated: June 10, 2003

Received: June 12, 2003

Dear Mr. Glover:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-46. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Susan Runner, DDS, MA

Interim Director

Division of Anesthesiology, General Hospital, Infector Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

ENCLOSURE B

Indications for Use

Page 1 of 1

<u>KO31814</u> (To be assigned)

Device Name:

StatCO₂meter™

Intended Use:

The Mercury Medical StatCO₂meter™ End Tidal CO₂ Detector with Integrated Airway Pressure Manometer is intended to provide a semi-quantitative visualization of the CO₂ in the patient airway while measuring Airway and PEEP pressure. It is an adjunct in patient assessment, to be used in conjunction with other methods to determine clinical signs and symptoms by or on the order of a physician.

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Anesthesiology, General Hospital,

Infection Control, Dental Devices

510(k) Number: <u>K 631814</u>

Prescription Use _X_ (Per CFR 801.109)

or

Over-the-counter use ___