

OCT 22 2003

SECTION 8

SUMMARY OF SAFETY AND EFFECTIVENESS

**510(k) Summary of
Safety and Effectiveness**

Information supporting claims of substantial equivalence, as defined under the Federal Food, Drug and Cosmetic Act, respecting safety and effectiveness is summarized below. For the convenience of the Reviewer, this summary is formatted in accordance with the Agency's final rule "...510(k) Summaries and 510(k) Statements..." (21 CFR 807) and can be used to provide a substantial equivalence summary to anyone requesting it from the Agency.

MODIFIED DEVICE NAME: Modified VICRYL* Plus
Antibacterial (Polyglactin 910) Suture

PREDICATE DEVICES NAME: Coated VICRYL* Plus
Antibacterial (Polyglactin 910) Suture.

Device Description

Modified Coated VICRYL* Plus Antibacterial (Polyglactin 910) suture is a synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide. Modified Coated VICRYL* Plus suture is coated with a mixture composed of equal parts of a copolymer of glycolide and lactide (Polyglactin 370) and calcium stearate and a small amount of an antimicrobial agent, Irgacare MP (*triclosan*).

Intended Use

Modified Coated VICRYL* Plus Antibacterial suture is intended for use in general soft tissue approximation and/or ligation, except for use in ophthalmic, cardiovascular and neurological tissues.

Indications Statement

Modified Coated VICRYL* Plus Antibacterial suture is indicated for use in general soft tissue approximation and/or ligation, except for ophthalmic, cardiovascular and neurological tissues.

* Trademark

**Technological
Characteristics**

The modified device has similar technological characteristics as the predicate devices. Like the currently marketed Coated VICRYL* Plus suture device, it is a sterile, braided synthetic absorbable suture that conforms to the USP Monograph for absorbable surgical sutures, except for diameter. Like the currently marketed Coated VICRYL Plus Antibacterial suture, the modified device contains Irgacare** MP, an antibacterial agent.

Performance Data

Non-clinical laboratory testing was performed demonstrating that the device conformed to the USP Monograph for absorbable surgical sutures. Additionally, in-vivo/in-vitro testing was provided showing that the device performed as intended and as claimed.

Conclusions

Based on the 510(k) summaries and 510(k) statements (21 CFR 807) and the information provided herein, we conclude that the modified device is substantially equivalent to the Predicate Devices under the Federal Food, Drug, and Cosmetic Act.

Contact

Rey Librojo
Senior Project Manager, Regulatory Affairs
ETHICON Products
ETHICON, Inc.
Rt. #22, West
Somerville, NJ 08876-0151

Date

August 1, 2003



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

OCT 22 2003

Mr. Rey Librojo
Regulatory Affairs
Ethicon, Inc.
Route 22 West
Somerville, New Jersey 08876

Re: K032420

Trade/Device Name: Modified Coated VICRYL* Plus Antibacterial (Polyglactin 910)
Synthetic Absorbable Suture

Regulation Number: 21 CFR 878.4493

Regulation Name: Absorbable poly(glycolide/L-lactide) surgical suture

Regulatory Class: II

Product Code: GAM

Dated: August 1, 2003

Received: August 7, 2003

Dear Mr. Librojo:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Mr. Rey Librojo

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", is written over the typed name.

Celia M. Witten, Ph.D., M.D.
Director

Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if known):

K032420

Device Name:

Modified Coated VICRYL* Plus Antibacterial (Polyglactin 910)
Synthetic Absorbable Suture

Indications for Use:

Modified Coated VICRYL* Plus Antibacterial suture is indicated for use in general soft tissue approximation and/or ligation, except for ophthalmic, cardiovascular and neurological tissues.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON
ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The Counter Use

(Optional Format 1-2-9G)

for Mark A. Miller
Division Sign-Off
Division of General Restorative
and Neurological Devices

510(k) Number K032420