

### Summary of Safety and Effectiveness information Special 510(k) – AEQUALIS Shoulder Fracture System K032679

**Regulatory authority:** Safe Medical Devices Act of 1990, 21 CFR 807.92

**1) Device name**

**Trade name:** *AEQUALIS Shoulder Fracture System*  
**Common name:** Total-Shoulder System and Hemi-Shoulder System  
**Classification name:** Shoulder joint metal/polymer semi-constrained cemented prosthesis

**2) Submitter**

Tornier S.A.  
B.P. 11 - Rue Doyen Gosse  
38330 Saint Ismier - France

**3) Company contact**

Tornier S.A.  
Mrs Mireille Lémery  
Regulatory affairs & Quality Engineer  
ZIRST - 161, rue Lavoisier  
38330 Montbonnot - France  
Tel: 00 33 4 76 61 38 98  
Fax: 00 33 4 76 61 35 33  
e-mail : [mireille.lemery@tornier.fr](mailto:mireille.lemery@tornier.fr)

**4) Classification**

**Device class:** Class II  
**Classification panel:** Orthopedic  
**Product code:** KWS  
§ 888.3660 Shoulder joint metal/polymer semi-constrained cemented prosthesis.

**5) Equivalent / Predicate device**

AEQUALIS Shoulder system, TORNIER SA (K952928)  
AEQUALIS Shoulder Fracture system, TORNIER SA (K994392, K003728)  
Modular Shoulder System, Wright Medical Technology Inc (K002683)



# TORNIER

## Implants Chirurgicaux

K032679

PAGE 2 of 2

### 6) Device description

The usual goal of total shoulder replacement and hemi-arthroplasty of the shoulder is to restore the shoulder joint to its best working condition and to reduce or eliminate pain. The *Aequalis Shoulder Fracture System* is intended to accomplish these goals. With the *Aequalis Shoulder Fracture System*, the natural glenoid elements of the shoulder may be conserved or replaced as warranted by the state of disease or injury. Thus the *Aequalis Shoulder Fracture System* is intended for use as a total shoulder replacement system, or as a hemi-shoulder. The modular nature of the system allows for the later conversion of a primary hemi-arthroplasty to a total shoulder replacement.

The present Device Modification submission consists in the addition of a long stem to each diameter of the previous range.

### 7) Materials

The stem is made of Titanium alloy (6Al-4V-Ti) according to ISO 5832-3. It is grit-blasted on its proximal part. The humeral head is made of Cobalt-Chromium alloy according to ISO 5832-7 or ISO 5832-12.

### 8) Indications

Traumatic or pathologic conditions of the shoulder resulting in fracture of the glenohumeral joint. Including humeral head fracture and displaced 3-or 4-part proximal humeral fractures.

Page 2/ page 2



TORNIER S.A  
161, rue Lavoisier  
38330 MONTBONNOT  
FRANCE

Tél. : 33 (0)4 76 61 35 00  
Fax : 33 (0)4 76 61 35 33

S.A. au capital de 1 800 000 F  
SIRET : 070 501 275 000 13  
R.C.S. : B 070 501 275  
CODE APE : 331 B

SIÈGE SOCIAL : B.P. 11 - rue du Doyen Gosse - 38330 SAINT-ISMIER - FRANCE



DEC 15 2003

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Mireille Lémery  
Regulatory Affairs and Quality Engineer  
Tornier S.A.  
ZIRST - 161, rue Lavoisier  
38330 Montbonnot  
FRANCE

Re: K032679  
Trade/Device Name: AEQUALIS Shoulder Fracture System  
Regulation Number: 21 CFR 888.3660  
Regulation Name: Shoulder joint metal/polymer semi-constrained cemented prosthesis  
Regulatory Class: II  
Product Code: KWS  
Dated: November 14, 2003  
Received: November 17, 2003

Dear Ms. Lémery:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", is written over the typed name.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Page 1 of 1

510(k) Number (if known): **K032679**

Device name: ***AEQUALIS Shoulder Fracture System***


**Indication for use:**

Traumatic or pathologic conditions of the shoulder resulting in fracture of the glenohumeral joint. Including humeral head fracture and displaced 3-or 4-part proximal humeral fractures.

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NECESSARY

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
Division Sign-Off)

Division of General, Restorative  
and Neurological Devices

510(k) Number K032679

Prescription use ✓

OR

Over-The-Counter Use \_\_\_\_\_

(Per 21 CFR 801.109)

(Optional format 1-2-96)