

FEB 23 2004

K033724

1.4 510(k) Summary of Safety and Effectiveness

Submitted by: Elizabeth J. Mason, Sr. Regulatory Affairs Specialist

Address: Nobel Biocare USA Inc.
22715 Savi Ranch Parkway
Yorba Linda, CA 92887

Telephone: (714) 282-4800, ext. 7830

Facsimile: (714) 998-9348

Date of Submission: November 21, 2003

Classification Name: Endosseous Implant (21 CFR 872.3640)

Trade or Proprietary
or Model Name: Nobel Biocare Centric Post

Legally Marketed Device(s): Replace® Scalloped Margin Implant System (K021584)
Nobel Perfect Implant System (K030257)

Device Description:

Nobel Biocare's Centric Post is a hollow, accessory component designed to enhance the performance of the implant/abutment system by ensuring proper alignment between an abutment and an endosseous implant, while additionally acting as a seal that prevents fluids from seeping into the implant interior.

The Nobel Biocare Centric Post does not contact any mucous membranes in the oral cavity because it rests inside the implant/abutment interface. The Centric Post can only be used as part of the implant system, not on its own. It is designed to remain inside the implant/abutment interface for as long as the implant remains in the patient's mouth.

These attributes of the Nobel Biocare Centric Post enhance the performance of the implant/abutment system, such as the Nobel Perfect Implant System, and make the component desirable for use in fabricating both final restorations and temporary restorations.

Indications for Use:

The Nobel Biocare Centric Post is a support component indicated for use as an alignment post to center an abutment and as a seal to prohibit fluids from seeping into the implant interior.

The Nobel Biocare Centric Post is a component within the implant system, which is intended for use in restoring the chewing function of fully edentulous and/or partially edentulous patients.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

FEB 23 2004

Ms. Elizabeth J. Mason
Senior Regulatory Affairs Specialist
Nobel Biocare USA, Incorporated
22715 Savi Ranch Parkway
Yorba Linda, California 92887

Re: K033724
Trade/Device Name: Nobel Biocare Centric Post
Regulation Number: 21 CFR 872.3640
Regulation Name: Endosseous Implant
Regulatory Class: III
Product Code: NHA
Dated: November 21, 2003
Received: November 26, 2003

Dear Ms. Mason:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



Chiu Lin, Ph., D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K033724

Device Name: Nobel Biocare Centric Post

Indications For Use:

The Nobel Biocare Centric Post is a support component indicated for use as an alignment post to center a temporary abutment and act as a seal to prohibit fluids from seeping into the implant interior during temporization.

The Nobel Biocare Centric Post is a component within the implant system, which is intended for use in restoring the chewing function of fully edentulous and/or partially edentulous patients.



(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K033724

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)