

Summary of Safety and Effectiveness Information [510(k) Summary] 3.0

SPONSOR:

Synthes (USA)

1690 Russell Road Paoli, PA 19301 (610) 647-9700

Contact: Lisa M. Boyle

DEVICE NAME:

Synthes (USA) Modification to Sternal Reconstruction System

CLASSIFICATION:

Class II § 21 CFR 888.3010: Cerclage, fixation

PREDICATE DEVICE:

Synthes Sternal Reconstruction System

DEVICE DESCRIPTION:

The Synthes modified Sternal Reconstruction System contains a 1.0 mm

stainless steel cable with detachable needle with ferrule The cable is 750

mm in length.

INTENDED USE:

Synthes (USA) Sternal Reconstruction System is intended for use in sternal

Comparative information presented supports substantial equivalence.

repair and reconstruction.

SUBSTANTIAL

EQUIVALENCE:





MAR - 2 2004

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Ms. Lisa M. Boyle Regulatory Associate Synthes (USA) 1690 Russell Road Paoli, Pennsylvania 19301

Re: K033816

Trade/Device Name: Modification to Synthes (USA) Sternal Reconstruction System

Regulation Number: 21 CFR 888.3010
Regulation Name: Bone fixation cerclage

Regulatory Class: II Product Code: JDQ Dated: December 8, 2003 Received: December 9, 2003

Dear Ms. Boyle:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Celia M. Witten, Ph.D., M.D.

Miriam C. Provost

Director

Division of General, Restorative and Neurological Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

2.0 Indications for Use Statement

		Page _1 of1
510(k) Number (if known): <u>Kc338</u>	16	
Device Name: Modification to Synthe	es (USA) Sternal Recons	truction System
Indications:		
Synthes (USA) Sternal Reconstruction System is intended for use in sternal repair and reconstruction.		
(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)		
Concurrence of CDRH, Office of Device Evaluation (ODE)		
Prescription Use (Per 21 CFR 801.109)	OR	Over-The-Counter Use_
Muram C. Provost (Division Sign-Off) Division of General, Restorative,		
and Neurological Devices		
510(k) Number.	K633816	