

APR 26 2004

December 24, 2003

K034063

Page 1 of 3

**510(k) SUMMARY AS REQUIRED BY SECTION 807.92(c)**

**Submitted by:**

Irvine Scientific Sales Co., Inc.  
2511 Daimler Street  
Santa Ana, CA 92705-5588

Telephone: (800) 437-5706  
Facsimile: (949) 261-6522

Contact: Wendell Lee

Date Submitted: November 21, 2003

**Device Identification:**

Trade Name:	MultiBlast Medium
Common Name:	In vitro embryo culture medium
Classification Name:	Reproductive Media (21 CFR, 886.6180)

**Predicate Device:**

Notice of Final Rule, 63 FR 48428, Docket number 97N-0335

**Description:**

MultiBlast Medium is a synthetic, defined medium, composed of a balance mixture of salts, amino acids, minerals and other nutrient substances designed to support embryonic growth and blastocyst development in vitro.

**Intended Use:**

MultiBlast Medium is intended for use as the second stage of a sequential in vitro embryo culture protocol. MultiBlast Medium has been developed to support the growth of human embryos from day three through day five, post-fertilization, including blastocyst formation.

**Technological Characteristics:**

After allowing the fertilized zygote to develop in vitro in a less complex, glucose- and phosphate-free culture medium (usually through day three, post-fertilization), the embryo is removed from the culture dish. It is placed into a fresh dish containing MultiBlast Medium, and protein supplementation. The dish is then returned to the incubator, and allowed to continue development, in vitro, until the desired stage of development has been achieved (usually day five post-fertilization). At that time, the embryo is removed from the medium, placed into a HEPES-buffered transport medium, and implanted into the patient.

**Performance Data:**

MultiBlast Medium is assayed by mouse embryo assay prior to release to market. This assay assures that the product is both functional for its intended use, the support of embryonic growth, and that no toxic components are present in the formulation. MultiBlast Medium has been used in a variety of clinical settings, for its intended use, for a number of years. In that time, the product has become the one of the standard media used as the second, more complex stage of a sequential media protocol.

**Additional Information:**

Mouse embryo testing will be performed as a condition of release for this product, as well as endotoxin and sterility testing. Results of all release assays performed will be reported on a lot-specific certificate of analysis, and will be indicated on the labeling.

**Conclusion:**

The conclusion from performance testing, as well as a review of the historical information contained in professional literature shows that MultiBlast Medium is suitable for its intended use, and meets the criteria outlined in the Notice of Final Rule, 63 FR 48428, Docket number 97N-0335.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

APR 26 2004

Wendell Lee, Pharm. D.  
Vice President  
Regulatory Affairs/Quality Systems  
Irvine Scientific Sales Co., Inc.  
2511 Daimler Street  
SANTA ANA CA 92705-5588

Re: K034063

Trade/Device Name: MultiBlast Medium  
Regulation Number: 21 CFR §884.6180  
Regulation Name: Reproductive media and supplements  
Regulatory Class: II  
Product Code: 85 MQL  
Dated: April 5, 2004  
Received: April 9, 2004

Dear Dr. Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of the letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97) you may obtain. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K034063

Device Name: MultiBlast Medium

Indications For Use: MultiBlast Medium is intended for use in the culture of human embryos to the blastocyst stage of development.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)


~~AND~~/OR

Over-The-Counter Use             
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number   K034063  

Page 1 of   1