

510(K) SUMMARY

APR - 7 2004

K040105(P.1 of 2)

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA and 21 CFR §807.92

1. Submitter's Name: YUE PFONG International Industrial Co., Ltd.
No.18, 23 Road, Taichung Ind. Zone, Taiwan
Address: (04)2359-0970
Phone: (04)2359-5876
Fax: Mr. Chang Yao-Chamr
Contact:
2. Device Name :
Trade Name: **YUE PFONG® Safety Syringe (3ml)**
Common Name: Safety Syringe (provided with needle)
Classification name Anti-Stick Syringe
3. Classification: Class II
Regulatory Number: 880.5860
Product Code: MEG
4. Predicate Device: ♦ TMD™ Safety Syringe (FA12 Series 3ml)with 510K number K022278
Marketed by Taiject Medical Device Co., Ltd.
5. Device Description: The **YUE PFONG® Safety Syringe (3ml)** is sterile, single-use, disposable , Non-reusable, Manual , Retractable, 3ml Piston Syringe, provided with needle attached in place., which is used for injection of fluids into the body. (Needle Spec. including 1 1/2" or shorter, 18-28 gauge). The **YUE PFONG® Safety Syringe (3ml)** consist of the following major components.
① Plunger , ② Barrel , ③ Piston, ④ adaptor ,⑤ O- Ring ,
⑥ Needle & Needle Hub, ⑦ Needle Sheath.
6. Intended Use: The **YUE PFONG® Safety Syringe (3ml)** is designed as an anti-stick syringe to reduce the risk of sharps injuries and the potential for syringe reuse and is a single use, disposable and manual retractable safety syringe which is intended for injection of medical fluids into the body.

7. Performance Summary: In terms of Physical specification, Chemical specification, Biological specification & Sterilization Specification, the device conforms to applicable standards included ISO 7864, ISO 7886-1, ISO 10993 series, ISO 11607-1, ISO 11135, USP Pyrogenic standards & related standards---etc.

8. Conclusions:

The **YUE PFONG® Safety Syringe (3ml)** has the same intended use and similar technological characteristics as the **TMD™ Safety Syringe --FA12 Series 3ml** market by **Taiject Medical Device Co., Ltd.** Moreover, bench testing & simulated use study contained in this submission demonstrate that any differences in their technological characteristics do not raise any new questions of safety or effectiveness. Thus, the **YUE PFONG International Industrial Co., Ltd.** is substantially equivalent to the predicate devices.



APR - 7 2004

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Yue Pfong International Industrial Company Limited
C/O Ms. Jennifer Reich
3892 South America West Trail
Flagstaff, Arizona 86001

Re: K040105

Trade/Device Name: YUE PFONG® Safety Syringe (3mL)
Regulation Number: 880.5860
Regulation Name: Piston Syringe
Regulatory Class: II
Product Code: MEG
Dated: March 4, 2004
Received: March 15, 2004

Dear Ms. Reich:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-56. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,



for, Chiu Lin, Ph.D.
Director
Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K#040105

Device Name: **YUE PFONG® Safety Syringe (3ml)**
YUE PFONG International Industrial Co., Ltd.

Indications For Use:

The YUE PFONG® Safety Syringe (3ml) is designed as an anti-stick syringe to reduce the risk of sharps injuries and the potential for syringe reuse and is a single use, disposable and manual retractable safety syringe which is intended for injection of medical fluids into the body.


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number: K040105

Page 1 of 1