

510(k) SUMMARY

[ENSPLY

NAME & ADDRESS:

DENTSPLY International
World Headquarters
Susquehanna Commerce Center
221 West Philadelphia Street
York, PA 17405-0872
(717) 845-7511 (voice)
(717) 849-4762 (fax)

www.dentsply.com

P. J. Lehn Telefax (717) 849-4343

CONTACT:

P. Jeffery Lehn

DATE PREPARED:

February 11, 2004

TRADE OR PROPRIETARY NAME: IDEAL® 1 ORTHODONTIC BAND CEMENT

CLASSIFICATION NAME:

Bracket adhesive and tooth conditioner resin (872.3750)

PREDICATE DEVICES:

Ideal® 1 Orthodontic Bracket Adhesive K033703

DEVICE DESCRIPTION: The IDEAL® 1 ORTHODONTIC BAND CEMENT is a light-cure two-part system: an orthodontic cement and a primer. A tooth conditioner is available if desired for the highest stress bonding cases.

The cement is packaged in syringes. No mixing is necessary. The primer is packaged in unit dose, self-dispensing capsules.

IN FENDED USE: The IDEAL® 1 ORTHODONTIC BAND CEMENT is indicated for boarding of orthodontic bands to natural and artificial tooth surfaces.

TECHNOLOGICAL CHARACTERISTICS: All of the components found in the IDEAL® 1 OF THODONTIC BAND CEMENT have been used in legally marketed devices.

We believe that the prior use of the components of the IDEAL® 1 ORTHODONTIC BAND CEMENT in legally marketed devices, the performance data provided, and the similarity of the prodicate to the new device support the safety and effectiveness of the IDEAL® 1 ORTHODONTIC BAND CEMENT for the intended uses.



APR - 2 2004

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Dentsply International Mr. P. Jeffery Lehn Director of Corporate Compliance & Regulatory Affairs World Headquarters Susquehanna Commerce Center 221 West Philadelphia Street York, Pennsylvania 17405-0872

Re: K040375

Trade/Device Name: IDEAL®1 Orthodontic Band Cement

Regulation Number: 872.3750

Regulation Name: Bracket Adhesive Resin and Tooth Conditioner

Regulatory Class: II Product Code: DYH Dated: February 11, 2004

Received: February 17, 2004

Dear Mr. Lehn:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-56. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html

Sincerely yours,

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital, Infection Control and Dental Devices Office of Device Evaluation Center for Devices and

Radiological Health

Enclosure

INDICATIONS FOR USE STATEMENT

(As Required by 21 CFR 807.87(e)

510(K) Number (if known): KOY	0375	
De rice Name: IDEAL® 1 ORTH	ODONTIC BAND	CEMENT
Indications for Use:		
Bonding of orthodontic bands to na	atural and artificial to	ooth surfaces
Prescription Use X (Part 21 CFR 801 Subpart D)	AND/OR	Over-The-Counter Use(21 CFR 807 Subpart C)
(I LEASE DO NOT WRITE BELOW TH	IIS LINECONTIN	NUE ON ANOTHER PAGE IF NEEDED)
Concurrence of CD	RH, Office of Devic	ee Evaluation (ODE)
(Division Sign- Division of And Infection Contr	Off) esthesiology, General rol, Dental Devices	Hospital,
510(k) Numbe	. (