

MAR 17 2004

K040510 Page 1/2

510(K) SUMMARY OF SAFETY AND EFFECTIVENESS

Common/Usual Name: Topical Hemostat

Product Trade Name: D-Stat Dry 3x3 Hemostatic Pad

Classification Name: Unclassified
Product Code FRO

Manufacturer: Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, MN 55369

Establishment Registration: 2134812

Contact: Gregory W. Sachs
Director of Regulatory Affairs

Performance Standards: No performance standards have been developed under section 514 for this device.

Device Description:

The D-Stat Dry 3x3 Hemostatic Pad consists of a:
Lyophilized pad consisting of thrombin, sodium carboxymethylcellulose and calcium chloride

The D-Stat Dry 3x3 Hemostatic Pad achieves its principal intended action (hemostasis) in the same manner as the other devices in the D-Stat Dry Product Family; it is applied directly over the source of bleeding therefore creating a physical barrier to blood flow through the application of adjunctive manual compression. The lyophilized components (thrombin, CMC, and calcium chloride) establish an environment in which a natural blood clot can build and form a physical barrier to bleeding. The thrombin facilitates hemostasis by enhancing the surface-activated clotting cascade through enzymatic cleavage and conversion of fibrinogen to fibrin.

Intended Use:

The Vascular Solutions D-Stat Dry 3x3 Hemostatic Pad is intended to be a topical compression bandage for use under the direction of a healthcare professional and is indicated as a trauma dressing for temporary control of moderate to severely bleeding wounds and for the control of surface bleeding from vascular access sites and percutaneous catheters or tubes.

Summary of Non-Clinical Testing:

Testing conducted included assessments of the physical properties of the lyophilized pad and packaging and the ability of the lyophilized components to achieve its intended use

(clot blood). The results of this battery of tests confirmed the suitability of the D-Stat Dry 3x3 Hemostatic Pad for its intended use.

Summary of Clinical Testing:

No clinical evaluations of this product for this use have been conducted.

Predicate Devices:

The intended use of the D-Stat Dry 3x3 Hemostatic Pad is a combination of the intended use of the:

- Vascular Solutions Inc., D-Stat Dry Hemostatic Bandage
- On Site Gas Systems, Inc., QuikClot,
- Marine Polymer Technologies, Inc. RDH Bandage

Conclusions:

The D-Stat Dry 3x3 Hemostatic Pad is substantially equivalent to Vascular Solutions Inc., D-Stat Dry Hemostatic Bandage, the On Site Gas Systems QuikClot and the Marine Polymer Technologies, Inc. RDH Bandage. The testing performed confirms that the D-Stat Dry 3x3 Hemostatic Pad will perform as intended.



MAR 17 2004

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Gregory W. Sachs
Director of Regulatory Affairs
Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, Minnesota 55369

Re: K040510

Trade/Device Name: Vascular Solutions D-Stat —Dry™ 3x3 Hemostatic pad
Regulatory Class: Unclassified
Product Code: FRO
Dated: February 26, 2004
Received: February 27, 2004

Dear Mr. Sachs:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,


for Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use Statement

510(k) Number: K040510

Device Name: Vascular Solutions D-Stat—Dry™ 3x3 Hemostatic Pad

Indications for Use:

The D-Stat Dry 3X3 is applied topically and is indicated as a trauma dressing for temporary control of moderate to severely bleeding wounds and for the control of surface bleeding from vascular access sites and percutaneous catheters or tubes.

Prescription Use ☒
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Miriam C. Provost
(Division Sign-Off)

**Division of General, Restorative
and Neurological Devices**

510(k) Number K040510

Page 1 of _____