

JUL 23 2004

Life Spine, LLC
NEO™ Cervical Plating System
510(k) Summary of Safety and Effectiveness

K040844
page 1 of 1

1. General Information

Submitter Information: Life Spine, LLC
11925 East 65th Street, Suite 4
Indianapolis, IN 46236
Phone: 317-826-7006 Fax: 317-826-7007

Company Registration Number: 3004499989

Date Summary Prepared: March 31, 2004

Classification Name: KWQ 888.3060 – Spinal Intervertebral Body Fixation
Orthosis. Class II

Common or Usual Name: Spinal Fixation System

Device Trade Name: NEO™ Anterior Cervical Plate System

2. Predicate Devices

Interpore Cross International Anterior Cervical Plate System (K002592) and the Howmedica Osteonics Corporation Reflex Anterior Cervical Plate System (K031702).

3. Intended Use

The NEO Anterior Cervical Plate System is indicated for use in temporary stabilization of the anterior spine from C2 to T1 during the development of cervical spinal fusions in patients with:

1. Degenerative disc disease, DDD (as defined by neck pain of discogenic origin with degeneration of disc confirmed by patient history and radiographic studies);
2. Spondylolisthesis
3. Trauma (including fractures or dislocations);
4. Spinal cord stenosis;
5. Deformity or curvatures (i.e. kyphosis, lordosis or scoliosis);
6. Tumors;
7. Pseudarthrosis;
8. Failed previous fusions.

WARNING: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

4. Device Description

The NEO Anterior Cervical Plate System consists of various sizes of anterior cervical bone plates, screws and screw locking tabs. Components are available in a variety of sizes to fit patient anatomy. All components are manufactured from implant grade titanium alloy 6Al-4V ELI per ASTM F-136. The NEO ACPS components will be supplied clean and "NON-STERILE".

5. Mechanical Testing

The results of the static and fatigue strength testing for the NEO Anterior Cervical Plate System demonstrate comparable mechanical properties to the predicate devices listed.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUL 23 2004

Michael S. Butler
President
Life Spine, LLC.
11965 East 65th Street, Suite 4
Indianapolis, Indiana 46236

Re: K040844
Trade/Device Name: NEOTM Anterior Cervical Plate System
Regulation Number: 21 CFR 888.3060
Regulation Name: Spinal intervertebral body fixation orthosis
Regulatory Class: II
Product Code: KWQ
Dated: June 29, 2004
Received: June 30, 2004

Dear Mr. Butler:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", with a stylized flourish at the end.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

INDICATIONS FOR USE STATEMENT

Page 1 of 1

510(k) Number (if known): K040844

Device Name: NEO™ Anterior Cervical Plate System

Indications for Use:

The NEO™ Anterior Cervical Plate System is intended for anterior interbody screw fixation of the cervical spine. This system is indicated for use in temporary stabilization of the anterior spine from C2 to T1 during the development of a cervical spinal fusions in patients with degenerative disc disease (as defined by neck pain of discogenic origin with degeneration of disc confirmed by patient history and radiographic studies); spondylolisthesis; trauma (including fractures or dislocations); spinal cord stenosis; deformity or curvatures (i.e. kyphosis, lordosis or scoliosis); tumors; pseudarthrosis; and / or failed previous fusions.

WARNING: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic, or lumbar spine.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE AS NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use: X
(Per 21 CFR 801.109)

Over-the-Counter Use: _____
(Optional Format 1-2-96)

[Signature]
(Division Sign-Off)
Division of General, Restorative,
and Neurological Devices

510(k) Number K040844