

JUN 1 5 2004

510(k) Summary

K040943 112

This summary of safety and effectiveness is submitted in accordance with the requirements of SDMA 1990 and 21 CFR 807.92

Establishment Registration Number: 2021898

Address of Manufacturer: Medtronic Neurosurgery  
125 Cremona Drive  
Goleta CA, 93117  
Telephone: (805) 968-1546 ext. 1770  
Fax: (805) 968-9336

Contact Person: Jeffrey Henderson

Date: April 8, 2004

Trade or Proprietary Name: Medtronic PS Medical Strata® Valves and Handtools

Common usual or Classification Name: Central Nervous System Fluid Shunts and Components, 21 CFR 882.5550

Description:

After surgical implantation, the Strata® Valve provides a non-invasive method to address changing patient needs in the management of hydrocephalus. The unique valve design allows the physician to adjust the pressure/flow performance level by using a magnetic Adjustment Tool included with the Strata® Valve adjustment kit without the need for radiographic confirmation.

The Strata® Handtools, include the Locator Tool, Indicator Tool and Adjustment Tool. All three tools are required to set or change the pressure/performance level setting of the Strata® Valve.

Intended Use:

The Strata® Valve is a shunt component designed to provide continued Cerebrospinal Fluid (CSF) flow from the ventricles of the brain into the right atrium of the heart or the peritoneal cavity. The design allows the physician to non-invasively adjust valve pressure/flow performance level pre- and post-implantation by using an external magnetic Adjustment Tool without the need for radiographic confirmation.

Predicate Device Identification:

The Strata® Valve is substantially equivalent to the Medtronic PS Medical Strata® Valve based on the same intended use, materials, design and dimensions.

**510(k) Summary****Device Testing:**

The Strata® Valve is substantially equivalent to the Medtronic PS Medical Strata® Valve based on successful performance testing, biocompatibility studies, and clinical study data.

Comparison of the Indicator Tool pressure/performance level setting and corresponding x-ray at time of Discharge and during follow-up, yielded the following results:

	X-ray & Handtool MATCH	X-ray & Handtool DIFFERENT
Discharge Data	98.6%	1.4%
Adjustment Data	98.8%	1.2%
TOTAL	98.7%	1.3%



JUN 15 2004

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Jeffrey Henderson  
Vice President, Quality & Regulatory Affairs  
Medtronic Neurosurgery  
125 Cremona Drive  
Goleta, California 93117-5500

Re: K040943  
Trade/Device Name: PS Medical Strata<sup>®</sup> Valve  
Regulation Number: 21 CFR 882.5550  
Regulation Name: Central nervous system fluid shunt and components  
Regulatory Class: II  
Product Code: JXG  
Dated: April 9, 2004  
Received: April 20, 2004

Dear Mr. Henderson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Mr. Jeffrey Henderson

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

  
for Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

## 2. MODIFIED INDICATIONS FOR USE STATEMENT

---

### Indications for Use

510(k) Number (if known): K040943

Device Name: PS Medical Strata® Valve

#### Indications For Use:

"The Strata® Valve is a shunt component designed to provide continued Cerebrospinal Fluid (CSF) flow from the ventricles of the brain into the right atrium of the heart or the peritoneal cavity. The design allows the physician to non-invasively adjust valve pressure/flow performance level pre- and post-implantation by using an external magnetic Adjustment Tool without the need for radiographic confirmation."

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

Miriam C. Provost  
(Division Sign-Off)

**Division of General, Restorative,  
and Neurological Devices**

**510(k) Number** K040943