K041213

JUL - 7 2004

## 510(k) Summary

In accordance with section 513(I) of the SMDA and as defined in 21 CFR Part 807.3 final rule dated December 14, 1994, this summary is submitted by:

Tyco Healthcare/Kendall 15 Hampshire Street Mansfield, MA 02048 Date Prepared: May 6, 2004

#### 1. Contact Person

Wei Zhao Senior Specialist, Regulatory Affairs Tyco Healthcare/Kendall Telephone: (508) 261-8404 Fax: (508)261-8461

#### 2. Name of Medical Device

Classification Name: Common or Usual Name:

Urological catheter and accessories
Urinary Drainage Catheter

# 3. Identification of Legally Marketed Device

The proposed Tyco Healthcare/Kendall *DOVER® RED RUBBER ROBINSON* Catheter is substantially equivalent in intended use, function and mode of operation to the Tyco Healthcare/Kendall **DOVER®** *RED RUBBER ROBINSON* Catheter, which is legally marketed prior to May 28, 1976.

#### 4. Device Description

The Tyco Healthcare/Kendall **DOVER® RED RUBBER ROBINSON** Catheter is a sterile, single use, intermittent urinary drainage catheter made from latex using dipping technology.

### 5. Device Intended Use

The product is intended for intermittent catheterization to drain urine from the urinary bladder. The product is intended for use on patients who are not capable of voluntary urination.

## 6. Product Comparison

The proposed *DOVER® RED RUBBER ROBINSON* Catheter is substantially equivalent to the predicate device in that each product has the same intended use and same physical,

Ray/243
Phile 242

functional and performance characteristics. The proposed and predicate devices are made from the same materials and have the same design.

# 7. Nonclinical Testing

Biocompatibility testing of the proposed device has demonstrated that it meets the requirements of guidelines presented in the 10993 ISO Standard, Part 1, with the FDA modified matrix presented in memorandum G95-1.

**End of Document** 



Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

# JUL - 7 2004

Ms. Wei Zhao Senior Specialist, Regulatory Affairs Tyco Healthcare The Kendall Company 15 Hampshire Street MANSFIELD MA 02048

Re: K041243

Trade/Device Name: Tyco Healthcare/Kendall DOVER® RED RUBBER ROBINSON Catheter

Regulation Number: 21 CFR §876.5130

Regulation Name: Urological catheter and accessories

Regulatory Class: II Product Code: 78 KOD Dated: June 22, 2004 Received: June 23, 2004

Dear Ms. Zhao:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of the letter:

| 8xx.1xxx                         | (301) 594-4591 |
|----------------------------------|----------------|
| 876.2xxx, 3xxx, 4xxx, 5xxx       | (301) 594-4616 |
| 884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx | (301) 594-4616 |
| 892.2xxx, 3xxx, 4xxx, 5xxx       | (301) 594-4654 |
| Other                            | (301) 594-4692 |
| Other                            | ` '            |

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97) you may obtain. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/dsma/dsmamain.html.

Sincerely yours,

Mancy C. Brogdon
Nancy C. Brogdon

Director, Division of Reproductive, Abdominal and Radiological Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

# Indication for Use Statement

510(k) Number (if known): K041243

#### Device Name:

Tyco Healthcare/Kendall DOVER® RED RUBBER ROBINSON Catheter

## Indications for Use:

The proposed device is intended for use in the drainage of urine from the urinary bladder. The product is intended for intermittent catheterization on patients who are not capable of voluntary urination. This product is not designed for use as an indwelling catheter.

Please DO NOT Write Below This Line - Continue On Another Page if Needed

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use \_X\_\_\_ OR OVER-The -Counter Use \_\_\_\_ (Per 21 CFR 801.109)

(Division Sign Off)

Division of Reproductive, Abdominal,

and Radiological Devices

510(k) Number 1942