

MAY 28 2004

K041283

1.4 510(k) Summary of Safety and Effectiveness

Submitted by: Elizabeth J. Mason
Sr. Regulatory Affairs Specialist

Address: Nobel Biocare USA Inc.
22715 Savi Ranch Parkway
Yorba Linda, CA 92887

Telephone: (714) 282-4800, ext. 7830

Facsimile: (714) 998-9348

Date of Submission: May 12, 2004

Classification Name: Porcelain Tooth (21 CFR 872.3920)

Trade or Proprietary
or Model Name: Procera® Bridge Zirconia

Legally Marketed Device(s): Procera® Copings & Pontic (K032562)

Device Description:

Nobel Biocare's Procera® Bridge Zirconia is a prefabricated device intended for use as the core structure of an artificial prosthesis for placement in the oral cavity in order to restore chewing function.

The Procera® Bridge Zirconia may be two, three, or four units and is precision milled. The Procera® Bridge Zirconia can be cemented or bonded to either natural or artificial tooth abutments. It is personalized according to the specific dimensions of the patient's abutments so the bridge precisely fits, and properly functions, in the patient's jaw.

Nobel Biocare's Procera® Bridge Zirconia is manufactured from one solid piece of densely sintered zirconia.

Indications for Use:

Nobel Biocare's Procera® Bridge Zirconia is indicated for use as the core structure of an artificial prosthesis for partially edentulous patients in the need of prosthetic oral reconstruction in order to restore chewing function. The Procera® Bridge Zirconia may be two, three, or four units and are cemented to natural or artificial tooth abutments.

1.5 Performance Standards

No special guidance documents, relative to this device, were found.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAY 28 2004

Nobel Biocare AB
C/O Ms. Elizabeth J. Mason
Senior Regulatory Affairs Specialist
Nobel Biocare USA, Incorporated
22715 Savi Sanch Parkway
Yorba Linda, California 92887

Re: K041283
Trade/Device Name: Procera[®] Bridge Zirconia
Regulation Number: 872.3920
Regulation Name: Porcelain Tooth
Regulatory Class: II
Product Code: ELL
Dated: May 12, 2004
Received: May 13, 2004

Dear Ms. Mason:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4613. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Chiu Lin', with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

1.3

Indications for Use

510(k) Number (if known):

K 041283

Device Name: Procera® Bridge Zirconia

Indications For Use:

Nobel Biocare's Procera® Bridge Zirconia is indicated for use as the core structure of an artificial prosthesis for partially edentulous patients in the need of prosthetic oral reconstruction in order to restore chewing function. The Procera® Bridge Zirconia may be two, three, or four units and are cemented to natural or artificial tooth abutments.


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Infection Control, Dental Devices

510(k) Number:

K 041283

Page 1 of 1

000007