

JUL 22 2004

Page 1/3

PROPRIETARY INFORMATION – LINVATEC CORPORATION

June 14, 2004

SUMMARY OF SAFETY AND EFFECTIVENESS

In accordance with the requirements of the Safe Medical Device Act of 1990 and 21 CFR 807.92, Linvatec Corporation is hereby submitting the Special 510(k) Summary of Safety and Effectiveness for the Super Revo® Herculine™ Suture Anchor 510(k) Number K041713

A. Submitter

Linvatec Corporation
11311 Concept Boulevard
Largo, Florida 33773-4908
Registration Number: 1017294

B. Company Contact

Elizabeth M. Paul
Manager, Regulatory Affairs
(727) 399-5234 Telephone
(727) 399-5264 FAX

C. Device Name

Trade Name: Super Revo® Herculine™ Suture Anchor

Common Name: Suture Anchor

Classification Names: Fastener, Fixation, Nondegradable, Soft Tissue, 21 CFR 888.3040

Proposed Class/Device: Class II

Product Code: MBI

PROPRIETARY INFORMATION – LINVATEC CORPORATION

Summary of Safety and Effectiveness
SuperRevo® Herculine™ Suture Anchor
510(k) # K041713
June 14, 2004

D. Predicate/Legally Marketed Devices

Super Revo® Suture Anchor
Linvatec Corporation

510(k) # K003984

E. Device Description

The Super Revo® Herculine™ Suture Anchor is a titanium suture anchor implant pre-threaded with undyed (white) and colored (white with blue stripe), nonabsorbable, braided, ultra-high molecular weight, polyethylene suture. It is provided preloaded onto a disposable driver with a stainless steel shaft and ABS handle. The Super Revo® Herculine™ Suture Anchor is supplied sterile and single use.

The modification described in this Special 510(K) is to supply the Super Revo® suture anchor with the ultra-high molecular weight polyethylene suture.

This modification does not affect the device's intended use, fundamental scientific technology or performance specifications.

F. Intended Use

The Super Revo® Herculine™ Suture Anchor is intended to be used for rotator cuff repairs in the shoulder either arthroscopically or in a mini-open technique.

G. Substantial Equivalence

The Super Revo® Herculine™ Suture Anchor is substantially equivalent in intended use, scientific technology and design to the Super Revo® Suture Anchor. Testing has been conducted to assure

PROPRIETARY INFORMATION – LINVATEC CORPORATION

that providing the suture anchor with the ultra-high molecular weight, braided, polyethylene suture does not raise any new issues regarding safety and effectiveness.

TeleFlex Medical, Fall River, MA supplies the suture which is approved for commercial distribution in the U.S. under 510(k) numbers K033654 (white) and K040472 (colored) .



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUL 22 2004

Ms. Elizabeth Paul
Manager, Regulatory Affairs
Linvatec Corp.
11311 Concept Boulevard
Largo, Florida 33773-4908

Re: K041713

Trade/Device Name: Super Revo[®] Herculine[™] Suture Anchor
Regulation Number: 21 CFR 888.3040, 21 CFR 878.5000
Regulation Name: Non-degradable soft tissue fastener; Polyethylene suture
Regulatory Class: II
Product Code: MBI, GAT
Dated: June 14, 2004
Received: June 23, 2004

Dear Ms. Paul:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Ms. Elizabeth Paul

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten", is written over the typed name.

Celia M. Witten, Ph.D., M.D.
Director

Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K041713

Device Name: Super Revo® Herculine™ Suture Anchor

Indications For Use: The Super Revo® Herculine™ Suture Anchor is intended to be used for rotator cuff repairs in the shoulder either arthroscopically or in a mini-open technique.


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of General, Restorative,
and Neurological Devices

Page 1 of 1

510(k) Number K041713