K 042691

NOV - 8 2005



Cook Incorporated

P.O. Box 489

Bloomington, IN 47402-0489 Phone: 800 468-1379 www.cookgroup.com

510(k) SUMMARY

Submitted By:

COOK INCORPORATED

750 Daniels Way, P.O. Box 489

Bloomington, IN 47402-0489

Contact Person:

Earl E. Knight III, MPA

Tel: (812) 339-2235 Fax: (812) 332-0281

Date Prepared:

September 28, 2004

Device:

Trade Name:

Vertefix[™] Vertebroplasty Procedure Set Filler, Bone Cement (For Vertebroplasty)

Common/Usual Name Proposed Classification:

Polymethylmethacrylate (PMMA) Bone Cement

21 CFR §888.3027, Class II

Product Code/Panel Code:

NDN—Orthopedic

Intended Use:

VertefixTM Radiopaque Bone Cement is indicated for the fixation of vertebral compression fractures during a vertebroplasty procedure. Painful vertebral compression fractures may result from osteoporosis, benign lesions (hemangioma), and malignant lesions (metastatic cancers, myeloma).

Predicate Devices:

Manufacturer	Device Name	510(k) #
Stryker Corp.	Spineplex TM Radiopaque Bone Cement	K032945
Kyphon, Inc.	KyphX [®] HV-R Bone Cement Model C01A	K041584

Device Description:

The VertefixTM Vertebroplasty Procedure Set contains the Vertebroplasty Injector Kit and VertefixTM Radiopaque Bone Cement. The bone cement consists of two separate, premeasured sterilized components: 20g polymer powder and 9.2g liquid monomer. The powder contains 30% barium sulfate as a radiopacifier.

Substantial Equivalence:

The subject device is similar with respect to intended use, chemical composition, and fundamental scientific technology of commercially available predicate devices in terms of section 510(k) substantial equivalence; any differences that may exist do not significantly affect the safety and effectiveness of the device.

Non-Clinical Test Data:

The results of these tests provide reasonable assurance that VertefixTM Radiopaque Bone Cement has been designed and tested to assure conformance to the requirements for its use as a bone cement. Testing conducted on the device includes: material and chemical analyses, mechanical testing, and biocompatibility.





NOV - 8 2005

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Thalia Brine Regulatory Affairs Specialist Cook Incorporated 750 Daniels Way, P.O. Box 489 Bloomington, Indiana 47402-1379

Re: K042691

Trade/Device Name: VertefixTM Vertebroplasty Procedure Set

Regulation Number: 21 CFR 888.3027

Regulation Name: Polymethylmethacrylate (PMMA) bone cement

Regulatory Class: II Product Code: NDN Dated: October 4, 2005 Received: October 6, 2005

Dear Ms. Brine:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Thalia Brine

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address http://www.fda.gov/cdrh/industry/support/index.html.

Sincerely yours,

Sar

Mark N. Melkerson

Acting Director

Division of General, Restorative and Neurological Devices
Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use Statement

510(k) Number (if know	(n): <u>KO426</u>	91		
Device Name: Vertefix TM Vertebroplasty Procedure Set				
Indications for Use:				
compression fractures du	ring a vertebroplasty pr osteoporosis, benign le	d for the fixation of vertebral occedure. Painful vertebral compression esions (hemangioma), and malignant		
Prescription Use <u>X</u> (Part 21 CFR 801 Subpar		Over-The-Counter Use(Part 21 CFR 807 Subpart C)		
PLEASE DO NOT WRITE	BELOW THIS LINE—CO	NTINUE ON ANOTHER PAGE IF NEEDED		
Concurren	ce of CDRH, Office of	Device Evaluation (ODE)		
-	7			
•	vision Sign-Off) ision of General, ¹	Restorative,		

and Neurological Devices

510(k) Number K 042691