



October 8, 2021

Vascular Solutions, Inc.  
Deborah Neymark  
VP, RA, Clinical Research & Quality Systems  
6464 Sycamore Court  
Minneapolis, Minnesota 55369

Re: K042937  
Trade/Device Name: Pronto Extraction Catheter  
Regulation Number: 21 CFR 870.5150  
Regulation Name: Embolectomy catheter  
Regulatory Class: Class II  
Product Code: QEZ

Dear Deborah Neymark:

The Food and Drug Administration (FDA) is sending this letter to notify you of an administrative change related to your previous substantial equivalence (SE) determination letter dated January 31, 2005. Specifically, FDA is updating this SE Letter as an administrative correction because FDA has created a new product code to better categorize your device technology.

Please note that the 510(k) submission was not re-reviewed. For questions regarding this letter please contact Gregory O'Connell, OHT2: Office of Cardiovascular Devices, (301) 796-6075, [Gregory.Oconnell@FDA.HHS.gov](mailto:Gregory.Oconnell@FDA.HHS.gov).

Sincerely,

Gregory W.  
O'connell -S

Digitally signed by Gregory W.  
O'connell -S  
Date: 2021.10.08 10:30:17 -04'00'

Gregory O'Connell  
Assistant Director  
DHT2C: Division of Coronary  
and Peripheral Intervention Devices  
OHT2: Office of Cardiovascular Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JAN 31 2005

Ms. Deborah Neymark  
Vice President, Regulatory Affairs  
Vascular Solutions, Inc.  
6464 Sycamore Court  
Minneapolis, MN 55369

Re: K042937  
Trade/Device Name: Pronto™ Extraction Catheter  
Regulation Number: 21 CFR 870.5150  
Regulation Name: Embolectomy Catheter  
Regulatory Class: II  
Product Code: DXE  
Dated: December 23, 2004  
Received: January 4, 2005

Dear Ms. Neymark:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set


Page 2 – Ms. Deborah Neymark

comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050. This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



 Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

### Indications for Use

510(k) Number (if known): K042937

Device Name: Vascular Solutions Pronto™ Extraction Catheter

#### Indications For Use:

The Vascular Solutions Pronto Extraction Catheter is indicated for the removal of fresh, soft emboli and thrombi from vessels of the arterial systems.

Prescription Use   X   AND/OR Over-The-Counter Use         
(Part 21 CFR 801 Subpart D) (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Donna D. Velmer  
Special Agent in Charge  
Office of Cardiovascular Devices  
Device K042937

Page 1 of   1  

3

K042937

## 510(K) SUMMARY OF SAFETY AND EFFECTIVENESS

JAN 31 2005

**Common/Usual Name:** Embolectomy Catheter

**Product Trade Name:** Pronto Extraction Catheter

**Classification Name:** Class II  
21 CFR 870.5150  
Product Code, DXE

**Manufacturer:** Vascular Solutions, Inc.  
6464 Sycamore Court  
Minneapolis, Minnesota 55369

**Establishment Registration:** 2134812

**Contact:** Deborah Neymark  
Vice President, Regulatory Affairs

**Performance Standards:** No performance standards have been developed under section 514 for this device.

### Device Description:

The Pronto extraction catheter is a dual lumen catheter with related accessories. The extraction lumen allows for the aspiration and removal of emboli/thrombi using the included syringe, extension line and stopcock. The catheter has a rounded distal tip with a protected, sloped opening of the extraction lumen to facilitate advancement of the catheter into the blood vessel and maximize extraction of emboli/thrombi through the extraction lumen. Incorporated within the catheter distal tip is a non-blood contacting radiopaque marker for fluoroscopic visualization. Two visible markers are located on the shaft of the catheter to provide guidance to the user concerning the relative position of the Pronto and the associated guide catheter. The catheter is a monorail design with a distal flexible region and a proximal stiff region. The catheter has an approximate outer diameter of 0.065 inches, allowing delivery through standard 6Fr. guide catheters. The smaller wire lumen of the catheter is able to accommodate guide wires that are  $\leq 0.014$ " in diameter. The catheter will be available in working lengths of 40 to 145 cm, in increments of 15 cm. The proximal end of the catheter incorporates a standard luer adapter to facilitate the attachment of the catheter to the included extension line, stopcock, and syringe. A 74 micron filter basket is included for assistance in filtering the blood removed during the procedure for laboratory analysis of thrombus.

### Intended Use:

The Pronto Extraction Catheter is indicated for the removal of fresh, soft emboli and thrombi from vessels in the arterial system.

**Summary of Non-Clinical Testing:**

Non-clinical testing of this product modification included assessments of the durability of the marker bands and the mid-shaft bond.

**Summary of Clinical Testing:**

No clinical evaluations of this product for this use have been conducted.

**Predicate Devices:**

Currently marketed Pronto Extraction Catheter (K032763)

**Conclusions:**

The modified Pronto is substantially equivalent to the currently marketed Vascular Solutions Inc. Pronto Extraction Catheter.