

JAN - 7 2005

3.0 510(k) Summary

Sponsor: Synthes (USA)
1690 Russell Road
Paoli, PA 19301
(610) 647-9700

Device Name: Synthes (USA) chronOS™

Classification: Class II, 21 CFR §888.3045
Filler, Calcium Sulfate Performed Pellets

Class II, 21 CFR §880.5860
Piston Syringe

Predicate Device: Synthes chronOS™
Orthovita, Inc. *Vitoss®* Scaffold Synthetic Cancellous Bone Void Filler and *Imbibe II* Syringe

Device Description: Synthes chronOS™ is a porous, osteoconductive, resorbable bone void filler made from β-Tricalcium Phosphate (TCP). Synthes chronOS™ is resorbed and replaced by new bone during the healing process. The pore structure of chronOS™ provides a matrix for the ingrowth of bone. chronOS™ is provided sterile in a double sterile pack. It is available in various forms including granules, blocks, wedges and cylinders.

Synthes chronOS™ may be packaged with a perfusion syringe that is used to mix the bone void filler with the patient's blood components and/or bone marrow.

Intended Use: **Synthes chronOS™** is indicated for use in bony voids or gaps that are not intrinsic to the stability of the bony structure. **chronOS™** is indicated for use in the treatment of bony defects created surgically or through traumatic injury.

Synthes chronOS™ is intended to be gently packed or placed into bony voids or gaps of the skeletal system (i.e. the extremities, spine, and pelvis) and may be combined with autogenous blood and/or bone marrow. Following placement into the bony void, **chronOS™** resorbs and is replaced with bone during the healing process.

Substantial Equivalence: Documentation is provided which demonstrates that Synthes chronOS™ is substantially equivalent* to other legally marketed Synthes devices.

*The term “substantial equivalence” as used in the 510(k) notification is limited to the definition of substantial equivalence found in the Food, Drug and Cosmetic Act, as amended and as applied under 21 CFR 807, Subpart E, under which a device can be marketed without pre-market approval or reclassification;. A determination of substantial equivalency under this notification is not intended to have any bearing whatsoever on the resolution of patent infringement suits or any other patent matters. No statements related to, or in support of substantial equivalence herein shall be construed as an admission against interest under the US Patent Laws or their application by the courts.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JAN - 7 2005

Ms. Kathy Anderson
Regulatory Manager
Synthes (USA)
1690 Russell Road
Paoli, Pennsylvania 19301

Re: K043045

Trade/Device Name: Synthes (USA) Chronos with Perfusion Syringe
Regulation Number: 21 CFR 888.3045, 21 CFR 880.5860
Regulation Name: Resorbable calcium salt bone void filler, Piston Syringe
Regulatory Class: II
Product Code: MQV, FMF
Dated: January 3, 2005
Received: January 4, 2005

Dear Ms. Anderson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Kathy Anderson

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Celia M. Witten', is positioned above the printed name.

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health



Indications for Use

510(k) Number (if known):

K043045

Device Name:

Synthes (USA) chronOS™

Indications:

Synthes chronOS™ is indicated for use in bony voids or gaps that are not intrinsic to the stability of the bony structure.

Synthes chronOS™ is indicated for use in the treatment of bony defects created surgically or through traumatic injury.

Synthes chronOS™ is intended to be gently packed or placed into bony voids or gaps of the skeletal system (i.e. the extremities, spine, and pelvis) and may be combined with autogenous blood and/or bone marrow. Following placement into the bony void, chronOS™ resorbs and is replaced with bone during the healing process.

Prescription Use ☒
(Per 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use ☐
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)

Division of General, Restorative,
and Neurological Devices

510(k) Number K043045

0004