

DEC 21 2004

K043280



510(k) Summary

Applicant/Sponsor: Walter Lorenz Surgical, Inc.
(A wholly owned subsidiary of Biomet, Inc.)
1520 Tradeport Drive
P.O. Box 18009
Jacksonville, FL 32229-8009
Establishment Registration Number: 1032347

Contact Person: Kacy Arnold, RN, MBA
Biomet Regulatory Specialist
Telephone: 574-267-6639
Fax: 574-372-1683

Proprietary Name: Mimix™ MP Bone Void Filler

Common Name: bone void filler

Classification Name: Methyl Methacrylate for Cranioplasty (882.5300)

Legally Marketed Devices

To Which Substantial

Equivalence Is Claimed: Mimix™ Bone Void Filler K990290

Device Description: The Mimix™ MP Bone Void Filler is packaged as separate, pre-measured powder and liquid components. The two components are designed to be mixed intraoperatively to produce a homogenous paste, which can then be applied to bone gaps or defects. Because of its apatitic nature, the material is highly biocompatible.

The powder component is a mixture of a ceramic calcium phosphate powder and sodium citrate dihydrate ($\text{Na}_3\text{C}_6\text{H}_5\text{O}_7 \cdot 2\text{H}_2\text{O}$). The liquid component is a solution comprised of anhydrous citric acid ($\text{C}_6\text{H}_8\text{O}_7$) and water (H_2O).

Indications for Use:

Mimix™ MP Bone Void Filler is a self-setting calcium phosphate cement indicated for the following craniofacial procedures

1. Repair of neurosurgical burr holes
2. Craniotomy cuts and other cranial defects
3. Augmentation or restoration of bony contour in the craniofacial skeleton area no larger than 25cm^2 .

Summary of Technologies: Mimix™ MP is substantially equivalent in materials and processing to the predicate device.

Non-Clinical Testing: Non-clinical testing demonstrated substantial equivalence between this device and the predicate device.

Clinical Testing: Clinical testing was not used to establish substantial equivalence.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

DEC 21 2004

Ms. Kacy Arnold, RN, MBA
Biomet Manufacturing Corp.
56 East Bell Drive
P.O. Box 587
Warsaw, Indiana 46582

Re: K043280

Trade/Device Name: Mimix™ MP Bone Void Filler
Regulation Number: 21 CFR 882.5300
Regulation Name: Methyl methacrylate for cranioplasty
Regulatory Class: II
Product Code: GXP
Dated: November 24, 2004
Received: November 26, 2004

Dear Ms. Arnold:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

for Miriam C. Provost

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K043280

Device Name: **Mimix™ MP Bone Void Filler**

Indications For Use:

Mimix™ MP Bone Void Filler is a self-setting calcium phosphate cement indicated for the following craniofacial procedures:

1. Repair of neurosurgical burr holes
2. Craniotomy cuts and other cranial defects
3. Augmentation or restoration of bony contour in the craniofacial skeleton area no larger than 25cm².

Miriam C. Provost
(Division Sign-Off)
**Division of General, Restorative,
and Neurological Devices**

510(k) Number K043280

Prescription Use XX
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use NO
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)