

MAR 18 2005

K043329 1 of 2

**510(k) Summary of Safety and Effectiveness for the  
Photo Therapeutics Limited Omnilux Revive and Blue Combination**

This 510(k) Summary of Safety and Effectiveness is being submitted in accordance with the requirements of the SMDA 1990 and 21 CFR 807.92.

**1. General Information**

Submitter:

Photo Therapeutics Limited  
Station House  
Stamford New Road  
Altrincham  
Cheshire WA14 1EP  
United Kingdom

Contact Person:

Steve Hutson  
Director of Engineering and Regulatory  
Affairs  
Photo Therapeutics Limited  
Station House  
Stamford New Road  
Altrincham  
Cheshire WA14 1EP  
United Kingdom

Summary Preparation Date:

November 17th, 2004

**2. Names**

Device Name:

Omnilux Revive and Blue Combination

Classification Name:

Laser Instrument, Surgical powered  
Product Code: GEX  
Panel: 79

**3. Predicate Devices**

Predicates for the Omnilux Revive and Omnilux Blue were detailed in K030426 and K030883 respectively, but the predicates are listed below for completeness.

Omnilux Revive is substantially equivalent to ...  
IPL Quantum SR manufactured by Lumenis, Inc. (K020839)  
Aurora SR manufactured by Syneron Medical Ltd.(K022266)  
EsteLux manufactured by Palomar Medical Technologies, Inc. (K020453)

00001

Omnilux Blue is substantially equivalent to ...  
Lumenis ClearLight (K013623).

#### **4. Device Description**

The Omnilux Revive and Blue is a combination of two sources of high spectral purity. They provide uniform or "hot-spot" free illumination. The two output are pre-tuned to one wavelength with a narrow spectral bandwidth. The output wavelength of Omnilux Revive is  $633 \pm 5$  nm, and Omnilux Blue is  $415 \pm 5$  nm. The Omnilux base unit contains the power supplies and the control unit. Attached to the base unit are three folding arms. The LED head can be attached to the end of the arms and then positioned for patient treatment. The control unit consists of an LCD and keyboard together with the control electronics. The user interface software allows the operator to access and control all device functions.

#### **5. Indications for Use**

The Omnilux Revive and Blue Combination is intended to emit energy in the red and blue region of the spectrum to treat dermatological conditions, specifically indicated to treat mild to moderate acne vulgaris.

#### **6. Performance Data**

Based upon an analysis of the overall performance characteristics for the device, Photo Therapeutics Limited believes that no significant differences exist between the 'Omnilux Revive and Blue Combination' and the previously approved Omnilux Revive (K030426) and Omnilux Blue (K030883) devices and their respective predicates. Therefore, the Omnilux Revive and Blue Combination raises no new issues of safety or effectiveness.

The new intended use is detailed in Appendix 1.



**DEPARTMENT OF HEALTH & HUMAN SERVICES**

Public Health Service

MAR 18 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Steve Hutson  
Director of Engineering and Regulatory Affairs  
Photo Therapeutics Ltd.  
Station House  
Stamford New Road  
Altrincham  
Cheshire WA14 1EP  
United Kingdom

Re: K043329  
Trade/Device Name: Omnilux Revive and Blue Combination  
Regulation Number: 21 CFR 878.4810  
Regulation Name: Laser surgical instrument for use in general and  
plastic surgery and in dermatology  
Regulatory Class: II  
Product Code: GEX  
Dated: February 1, 2005  
Received: February 7, 2005

Dear Mr. Hutson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Mr. Steve Hutson

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'M. Provost', with a stylized flourish at the end.

Miriam C. Provost, Ph.D.  
Acting Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

K043329

Device Name Omnilux Revive and Blue Combination

**Indications for Use:**

The Omnilux Revive and Blue Combination is intended to emit energy in the red and blue region of the spectrum to treat dermatological conditions, specifically indicated to treat mild to moderate acne vulgaris.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE  
IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over The Counter Use ☐



(Optional Format 1-2-96)

Director, Center for Devices and Radiological Services

K043329