

VII. 510(k) Summary

In accordance with the Safe Medical Devices Act (SMDA) of 1990 and Title of the Code of Federal Regulations Part 807 (21 CFR §807), and in particular §807.92, the following summary of safety and effectiveness information is provided:

A. Submitted by

Brian Burkinshaw
LDR SPINE USA
108 Wild Basin Rd.
Austin, Texas 78746
Telephone: (512) 437-4311
Date Prepared: November 16, 2004

B. Device Name

Trade or Proprietary Name: BF+
Common or Usual Name: Bone Void Filler
Classification Name: Unclassified

C. Predicate Devices

The subject device is substantially equivalent to similar previously cleared devices.

D. Device Description

The BF+ is a synthetic, resorbable calcium phosphate bone void filler. It is an osteoconductive material which provides a porous scaffold upon which bone formation can occur. The interconnected porosity ranges from 60 to 80% with a pore size range of 200 to 500µm. The device is available in a variety of shapes and sizes.

E. Intended Use

BF+ is indicated for filling bone voids or defects of the skeletal system (such as the extremities, spine and the pelvis) that are not intrinsic to the stability of the bony structure. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. BF+ is a bone graft substitute that resorbs and is replaced with bone during the healing process.

F. Substantial Equivalence

Data was provided which demonstrated the BF+ to be substantially equivalent to previously cleared devices. The substantial equivalence is based upon equivalence in indications for use, design, material, and function.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JAN - 7 2005

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Brian Burkinshaw
Dir. Technology Solutions
LDR Spine USA
108 Wild Basin Road
Austin, Texas 78746

Re: K043347
Trade/Device Name: BF+ Bone Substitute
Regulation Number: 21 CFR 888.3045
Regulation Name: Resorbable calcium salt bone void filler device
Regulatory Class: II
Product Code: MQV
Dated: November 30, 2004
Received: December 6, 2004

Dear Mr. Burkinshaw:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Mr. Brian Burkinshaw

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

for Miriam C. Provost

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

A. Indications for Use510(k) Number (if known): K043347Device Name: BF+**Indications for Use:**

BF+ is indicated for filling bone voids or defects of the skeletal system (such as the extremities, spine and the pelvis) that are not intrinsic to the stability of the bony structure. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. BF+ is a bone graft substitute that resorbs and is replaced with bone during the healing process.

Prescription Use ☒
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Miriam C. Provost
(Division Sign-Off)**Division of General, Restorative,
and Neurological Devices****510(k) Number** K043347