

510(K) SUMMARY *K043352 1/1*

MAR 31 2005

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA and 21 CFR §807.92

1.0 Submitter's Name: **Gocanning Technology Ltd.**

Address: 6F-1, No. 188 , Sec. 3, Tatung Rd., Hsichih City , Taipei , Taiwan
Phone: TEL: 001-886-2-86471618
Fax: FAX: 001-886-2-86472653
Contact: Mr. TSUNG-Kun Su, General Manager

2.0 Device Name

Trade Name: **GOCANNING GLS-700 e-Pulse Analysis Monitor**
Model No.: **GLS-700**

Common Name: **Non-Invasive Blood Pressure Monitor**

Classification name: System , Measurement , Blood-Pressure , Non-Invasive

3.0 Classification: Class II

4.0 Predicate Device: The predicate device is the Rossmax Automatic Blood Pressure Monitor(K021225) marketed by ROSSMAX INTERNATIONAL LTD.

5.0 Device Description: **GOCANNING GLS-700 e-Pulse Analysis Monitor** is designed to measure the systolic and diastolic blood pressure, and pulse rate(heart of an individual).

6.0 Intended Use: **GOCANNING GLS-700 e-Pulse Analysis Monitor** is intended to measure human's Systolic, Diastolic blood pressure and pulse rate. All values can be read out in one LCD DISPLAY.

7.0 Performance Summary: In terms of operating specification, Safety & EMC requirements, the device conforms to applicable standards included EN-1060-1, EN-1060-3, ANSI/AAMI SP-10, IEC 60601-1 and IEC 60601-1-2 requirements.

8. Conclusions:

The **GOCANNING GLS-700 e-Pulse Analysis Monitor** have the same intended use and similar technological characteristics as Rossmax Automatic Blood Pressure Monitor(K021225) marketed by ROSSMAX INTERNATIONAL LTD. Moreover, bench testing contained in this submission and clinical testing supplied demonstrate that any differences in their technological characteristics do not raise and new questions of safety or effectiveness. Thus, the **GOCANNING GLS-700 e-Pulse Analysis Monitor** is substantially equivalent to the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

MAR 31 2005

Gocanning Technology, Ltd.
c/o Ms. Jennifer Reich
Harvest Consulting Corp.
3892 South America Trail West
Flagstaff, AZ 86001

Re: K043352

Trade Name: GLS-700 e-Pulse Analysis Monitor
Regulation Number: 21 CFR 870.1130
Regulation Name: Noninvasive Blood Pressure Measurement System
Regulatory Class: II (two)
Product Code: DXN
Dated: March 10, 2005
Received: March 15, 2005

Dear Ms. Reich:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 -- Ms. Jennifer Reich

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,


Bram D. Zuckerman, M.D.
Director

Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K043352

Device Name: **GLS-700 e-Pulse Analysis Monitor**
Gocanning Technology Ltd.

Indications For Use:

GLS-700 e-Pulse Analysis Monitor is used to measure automatically human's Systolic, Diastolic blood pressure and pulse rate.

The intended for use of this over-the-counter device is for adult of age 18 and above.


Prescription Use _____
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use V
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Cardiovascular Devices
510(k) Number K043352