

FEB 16 2005

1C 093387

510(k) SUMMARY

**Submitter's
Information**

VasCon LLC
9344 NW 13 Street
Miami, Florida 33172 USA
Telephone: 1.305.477.2406
Fax: 1.305.592.6605
Contact: Stephen F. Vadas, Ph.D.

Preparation Date

08 December 2004

Name of Device

Common Name: Guiding Catheter
Classification Name: Percutaneous Catheter
Trade Name: VasCon Polaris™ Guiding Catheter

Predicate Device

VasCon Polaris™ Guiding Catheters, 510(k) K033633

Intended Use

The guiding catheter is intended for use for intravascular introduction of interventional / diagnostic devices into the coronary or peripheral vascular systems.

**Device
Description and
Summary of
Technological
Characteristics**

The VasCon Polaris™ Guiding Catheters are single lumen catheters, incorporating a nylon body reinforced with a stainless steel wire braid, a PTFE lubricious inner lumen, and a soft radiopaque tip to reduce potential vessel injury. They are available in 5Fr, 6Fr, 7Fr, and 8Fr, 100 cm in length, and in a variety of shapes. The technological characteristics are equivalent to the predicate device.

**Testing
Summary**

Laboratory testing has been performed on the VasCon Polaris™ Guiding Catheters to assure compliance to the specifications.

Conclusions

The testing as discussed above demonstrate that, like the predicate device, the VasCon Polaris™ Guiding Catheters are safe and effective for its intended use.

000016



FEB 16 2005

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Stephen F. Vadas, Ph.D., P.E.
Vice President, Product Assurance and Regulatory Affairs
VasCon LLC
9344 NW 13 Street
Miami, FL 33172

Re: K043387
Trade/Device Name: VasCon Polaris™ Guiding Catheter
Regulation Number: 21 CFR 870.1250
Regulation Name: Percutaneous catheter
Regulatory Class: II
Product Code: DQY
Dated: December 8, 2004
Received: December 9, 2004

Dear Dr. Vadas:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set


Page 2 – Stephen F. Vadas, Ph.D., P.E.


forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



 Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K043387

Device Name: VasCon Polaris™ Guiding Catheter

Indications For Use: The guiding catheter is intended for use for intravascular introduction of interventional / diagnostic devices into the coronary or peripheral vascular systems.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Denise R. Kachner
(Division 1 Sign-Off)
Division of Cardiovascular Devices

510(k) Number K043387