

**510(k) Summary****Premarket Notification, Section 510(k)****LDR Spine USA****JUNE 13, 2005**

Regulatory Authority: Safe Medical Devices Act of 1990, 21 CFR 807.92

**1. Device Name:**Trade Name: **MC+ Partial VBR****JUN 30 2005****Common**

Name(s): Vertebral Body Replacement

**Classification**

Name(s): Vertebral Body Replacement (MQP)

**2. Establishment Name & Registration Number:**Name: **LDR Spine USA**

Number: 3004903783

**3. Classification(s):**

Sec. 888.3060 Spinal intervertebral body fixation orthosis.

(a) **Identification.** A spinal intervertebral body fixation orthosis is a device intended to be implanted made of titanium. It consists of various vertebral plates that are punched into each of a series of vertebral bodies. An eye-type screw is inserted in a hole in the center of each of the plates. A braided cable is threaded through each eye-type screw. The cable is tightened with a tension device and it is fastened or crimped at each eye-type screw. The device is used to apply force to a series of vertebrae to correct "sway back," scoliosis (lateral curvature of the spine), or other conditions.

Device Class: Class II for the requested indications

Classification Panel: Orthopaedic and Rehabilitation Devices Panel

Product Code(s): MQP

**4. Equivalent Predicate Device:****LDR Spine USA** believes that the **MC+** is substantially equivalent to the following:

K042268 – EBI CAS Spine Spacer System

K043206 – Pioneer Vertebral Spacer

K032064 – CPOD / LPOD VBR System

K033109 – Ellys and Aurys VBR

Equivalence can be seen in the design, material composition, surgical technique and intended use.

**5. Device Description:**

The MC+ vertebral implants consist of a series of square and rectangular D-shaped implants. The device is used in pairs as a partial vertebral body replacement in the thoracolumbar spine (from T1 to L5). The device is offered in seven different configurations to better approximate the anatomical variation observed in different vertebral levels and/or patient anatomy.

The LDR Spine USA, MC+ is comprised of a variety of components fabricated and manufactured from PolyEtherEtherKetone (PEEK) as described by ASTM F-2026. This material is utilized due to its radiolucent properties, which aid the surgeon in determining if fusion in the operative site has occurred.

The superior and inferior surfaces of the construct have a pattern of teeth to provide increased stability and to help prevent movement of the device. The upper surface of the implant is arched in the anterior/posterior direction for increased stability. Tantalum wire markers (ASTM F-560) are inserted into components to give surgeons a visual aid in determining the location of the implant, both inter and post-operatively. The device comes with optional anchoring clips made of titanium alloy.

**Materials:** all implants are made from implant grade PolyEtherEtherKetone polymer (PEEK) with tantalum alloy position markers and titanium alloy as indicated in the table below:

PEEK Optima LT	USP Class VI ASTM F-2026	ISO 10993
Tantalum	ASTM F-560	ISO 5832-3
Titanium	ASTM F-136,92	ISO 5832-3

**Indications for Use.** The MC+ is indicated for use for partial replacement (i.e.: partial vertebrectomy) of a diseased vertebral body that has been resected or excised due to tumor or trauma/fracture in order to achieve anterior decompression of the spinal chord and neural tissues, and to restore the height of a collapsed vertebral body. The components are intended for use in pairs for use as a partial vertebral body replacement in the thoracolumbar spine (from T1 to L5) and are intended for use with supplemental internal fixation such as the Easyspine System to properly utilize this system.

**Testing Summary.** Fatigue and static testing is complete. Samples were tested according to accepted engineering and scientific principals. Test results demonstrate that the system can be expected to perform in a manner equivalent to the comparison device.

**6. Applicant Name & Address:**

**LDR Spine USA**  
4030 W. Braker Ln., Ste. 360  
Austin, TX 78759  
**Office:** (512) 344-3333  
**Fax:** (512) 344-3350

**7. Company Contact:**

Mr. Edward E. Newton  
Dir. Reg. and Clinical Affairs  
**LDR Spine USA**  
4030 W. Braker Ln., Ste. 360  
Austin, TX 78759  
**Office:** (512) 344-3316  
**Fax:** (512) 344-3350

**8. Submission Correspondent:**

Mr. Brian Burkinshaw  
Dir. Innovation & Technology Solutions  
**LDR Spine USA**  
4030 W. Braker Ln., Ste. 360  
Austin, TX 78759  
**Office:** (512) 344-3304  
**Fax:** (512) 344-3350

**9. Performance Standards:**

United States Food and Drug Administration mandated performance standards for this device do not exist. Various voluntary performance standards are utilized. Voluntary standards utilized include ASTM, LDR Spine, USA Standard Operating Procedures, vendor & process certification and

qualification procedures, Quality Systems Regulations, ISO materials standards and ISO 13485 series quality regulations.

**LDR Spine USA** also meets appropriate general controls authorized under Sections 501, 502, 510, 516, 518, 519, and 520 of the Food, Drug, and Cosmetic Act.

#### 10. **Storage, Packaging & Sterilization Information:**

The implantable portions of the **MC+** are supplied "**STERILE**". The sterilization process is radiation and the selected protocol has been validated. The minimum Sterility Assurance Level (SAL) of at least  $10^{-6}$ .

The instruments are supplied non-sterile and must be sterilized prior to use. The recommended sterilization process for the instruments is high temperature steam autoclave sterilization. The referenced sterilization cycle produces a Sterility Assurance Level (SAL) of at least  $10^{-6}$ .

**The validated cycle is:**

Method: Steam  
Cycle: Gravity  
Temperature: 270°F (134°C)  
Exposure Time: 18 minutes

All packages containing implants or instruments should be intact upon receipt. Damaged packaging may indicate the presence of unsafe product. If the package or product is damaged, the product should not be used and should be returned. Product must be handled, stored and opened in such a way that it is protected from inadvertent damage or contamination. When used, the instruments must be placed into use following cleaning, sterilization and accepted surgical sterile technique.

#### 11. **Summary Comparison Table:**

Feature Comparison Table:

FEATURE	MC+	CAS Spine Spacer System	Pioneer Vertebral Spacer	CPOD LPOD VBR System	Ellys and Aurys VBR	SE?
<b>Indications for Use:</b>	The MC+ is indicated for use for partial replacement (i.e.: partial vertebrectomy) of a diseased vertebral body that has been resected or excised due to tumor or trauma/fracture in order to achieve anterior decompression of the spinal chord and neural tissues, and to restore the height of a collapsed vertebral body. The components are intended for use in pairs for use as a partial vertebral body replacement in the thoracolumbar spine (from T1 to L5) and are intended for use with supplemental internal fixation such as the Easyspine System to properly utilize this system.	Same	Same	Same	Same	YES
<b>Design:</b>	Cage	Cage	Cage	Cage	Cage	YES
<b>Sterile:</b>	Yes	Yes	Yes	Yes	Yes	YES

<b>Material:</b>	PEEK w/ tantalum markers	Titanium	Peek	Titanium	PEEK w/ titanium markers	YES
<b>Instruments:</b>	Specialized instruments required. Reusable. Autoclave before use.	Same	Same	Same	Same	YES
<b>Bone Graft:</b>	Required	Same	Same	Same	Same	YES
<b>Secondary Fixation:</b>	Required	Same	Same	Same	Same	YES
<b>How used:</b>	Pairs	Singly	Singly	Pairs	Singly and Pairs	YES
<b>K Number:</b>	K043479	K042268	K043206	K032064	K033109	YES
<b>Manufacturer:</b>	LDR Spine USA	EBI	Pioneer Surgical Technology	Thcken Spine	Scient'x	YES
<b>Product Code:</b>	MQP	MQP	MQP	MQP	MQP	YES



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUN 3 0 2005

Mr. Brian Burkinshaw  
Director Technology Solutions  
LDR Spine USA Incorporated  
4030 West Braker Lane, Suite 360  
Austin, Texas 78759

Re: K043479  
Trade/Device Name: MC+ Partial Vertebral Body Replacement (PVBR)  
Regulation Number: 21 CFR 888.3060  
Regulation Name: Spinal intervertebral body fixation orthosis  
Regulatory Class: II  
Product Code: MQP  
Dated: June 12, 2005  
Received: June 16, 2005

Dear Mr. Burkinshaw:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Brian Burkinshaw

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Miriam C. Provost". To the left of the signature is a small, stylized mark that looks like "Jza".

Miriam C. Provost, Ph.D.  
Acting Director  
Division of General, Restorative,  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number: **K043479**Device Name(s): **MC+ Partial Vertebral Body Replacement (PVBR)****Indications for Use:**

The MC+ is indicated for use for partial replacement (i.e.: partial vertebrectomy) of a diseased vertebral body that has been resected or excised due to tumor or trauma/fracture in order to achieve anterior decompression of the spinal chord and neural tissues, and to restore the height of a collapsed vertebral body. The components are intended for use in pairs for use as a partial vertebral body replacement in the thoracolumbar spine (from T1 to L5) and are intended for use with supplemental internal fixation such as the Easyspine System to properly utilize this system.

Prescription Use   X  

OR

Over-The-Counter Use       

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NECESSARY

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)

Division of General, Restorative

(Per 21 CFR 801.109) and Neurological Devices

(Optional format 1-2-96)

510(k) Number K043479