

91-2

2010年10月
 第10期
 第10期

[illegible][illegible][illegible][illegible][illegible][illegible][illegible][illegible]

1443494 0352
leaving the catheter securely seated in the bone. The catheter has a standard Luer lock that permits attachment of standard syringes and IV tubing for administration of drugs and fluids.

Indications for Use:

The PD-IO Disposable Intraosseous Infusion Needle and Handle provides intraosseous access in the proximal tibia, as an alternative to IV access during emergencies. The device is for use in pediatric patients.

Summary of the technological characteristics of the VidaCare PD-IO Disposable Intraosseous Infusion Needle and Handle compared to the predicate devices:

The predicate Cook, Disposable Intraosseous Infusion Needles K913258, VidaCare Corporation EZ-IO K032885 and the WaisMed Ltd BIG Pediatric K022415 were compared in the following areas and found to have similar technological characteristics and to be equivalent to the PD-IO Disposable Intraosseous Infusion Needle.

- Indications for use
- Design features
- Needle design
- Technique
- Sterility
- Biocompatibility
- Anatomical site
- Where used
- Standards met
- Target population
- Mechanical safety



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

FEB 18 2005

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

VidaCare Corporation
C/O Mr. Greg Holland
Regulatory Specialist
Regulatory Specialists, Incorporated
3722 Avenue Sausalito
Irvine, California 92606

Re: K043490

Trade/Device Name: PD-IO Disposable Intraosseous Infusion Needle and Handle
Regulation Number: 21 CFR 880.5570
Regulation Name: Hypodermic Single Lumen Needle
Regulatory Class: II
Product Code: FMI
Dated: December 16, 2004
Received: December 17, 2004

Dear Mr. Holland:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Chiu Lin", with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,
Infection Control and Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Page 1 of 1

510(k) Number (if known): K043490

Device Name: PD-IO Disposable Intraosseous Infusion Needle and Handle

Indications For Use:

The PD-IO Disposable Intraosseous Infusion Needle and Handle provides intraosseous access in the proximal tibia, as an alternative to IV access during emergencies. The device is for use in pediatric patients.

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____
(Optional Format 1-2-96)

(PLEASE DO NOT WRITE BELOW THIS LINE—CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Chetan V. Mehta

(Division Sign-Off)
Division of Anesthesiology, General Hospital,
Injection Control, Dental Devices

510(k) Number: K043490