

JAN 31 2005

**510(k) Summary – C.f.a.s. (Calibrator for Automated Systems) PUC (Proteins in Urine/CSF); Precinorm ® PUC and Precipath ® PUC**

---

**Introduction** According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence

---

**Submitter name, address, contact** Roche Diagnostics  
9115 Hague Rd  
Indianapolis IN 46250  
(317) 521-3723

Contact person: Theresa M. Ambrose

Date prepared: January 5, 2005

---

**Device 1 Name** Proprietary name: Roche Diagnostics C.f.a.s. (Calibrator for automated systems) PUC (Proteins in Urine/CSF)

Common name: C.f.a.s. PUC

Classification name: Calibrator, multi-analyte mixture

---

**Device 2 Name** Proprietary name: Roche Diagnostics Precinorm ® PUC (Proteins in Urine/CSF) and Precipath ® PUC (Proteins in Urine/CSF)

Common name: Precinorm ® PUC / Precipath® PUC

Classification name: Multi-analyte controls, all kinds (assayed and unassayed)

---

*Continued on next page*

## 510(k) Summary, continued

<b>Device 1 description</b>	C.f.a.s. PAC is a liquid, ready-for-use calibrator consisting of a buffered aqueous solution with biological materials added as required to obtain desired component levels. Values for constituent analytes are provided in the product labeling.
<b>Device 2 description</b>	Precinorm ® PUC/ Precipath ® PUC is a liquid ready-for-use control based on a buffered aqueous solution. Concentrations of control components have been adjusted to represent normal and pathological ranges. Values for constituent analytes are provided in the product labeling.
<b>Device 1 Intended use</b>	C.f.a.s. PUC is for use in the calibration of quantitative Roche methods on Roche clinical chemistry analyzers as specified in the enclosed value sheet.
<b>Device 2 Intended use</b>	Precinorm ® PUC/ Precipath ® PUC is for use in quality control by monitoring accuracy and precision for the quantitative methods as specified in the enclosed value sheet
<b>Device 1 Substantial Equivalence</b>	For C.f.a.s. PUC, Roche claims substantial equivalence to the currently marketed currently marketed C.f.a.s. PUC, K040264.
<b>Substantial equivalence comparison – Device 1</b>	The table below compares C.f.a.s. PUC with its predicate device (currently marketed C.f.a.s. PUC, K040264).

*Continued on next page*

## 510(k) Summary, continued

Characteristic	C.f.a.s. PUC (Predicate device, K040264)	C.f.a.s. PUC (Modified Device)
Intended Use	C.f.a.s. PUC is for use in the calibration of quantitative Roche methods on Roche clinical chemistry analyzers as specified in the enclosed value sheet.	Same
Format	Liquid ready-for-use calibrator based on a buffered aqueous solution. The concentrations of the calibrator components have been adjusted to ensure optimal calibration of the appropriate Roche methods on clinical chemistry analyzers.	Same
Stability	<u>Unopened</u> Stable at 2-8°C until expiration date <u>Opened:</u> Stable at 2 to 8°C for 4 weeks	Same
Level	Single level	Same
Constituent Analytes with Assigned Values	<ul style="list-style-type: none"> <li>• Albumin</li> <li>• Total protein</li> </ul>	<ul style="list-style-type: none"> <li>• Albumin</li> <li>• Total protein</li> <li>• Immunoglobulin G</li> </ul>

### Device 2 Substantial Equivalence

For Precinorm ® PUC and Precipath ® PUC, Roche claims substantial equivalence to the currently marketed Roche Diagnostics Precinorm ® PUC and Precipath ® PUC (K041812).

### Substantial Equivalence – Device comparison

The table below compares Precinorm ® PUC / Precipath® PUC with the predicate device (currently marketed Precinorm ® PUC / Precipath® PUC).

*Continued on next page*

## 510(k) Summary, continued

Characteristic	Precinorm® PUC / Precipath® PUC (Predicate device, K040264)	Precinorm® PUC / Precipath® PUC (Modified Device)
Intended Use	For use in quality control by monitoring accuracy and precision for the quantitative methods as specified in the enclosed value sheet	Same
Format	Liquid ready-for-use control based on a buffered aqueous solution. Concentrations of control components have been adjusted to represent normal and pathological ranges.	Same
Stability	Unopened Stable at 2-8°C until expiration date Opened: Stable at 2 to 8°C for 4 weeks	Same
Constituent Analytes with Assigned Values	<u>Precinorm</u> <ul style="list-style-type: none"> <li>• Albumin</li> <li>• Creatinine</li> <li>• Total Protein</li> </ul> <u>Precipath</u> <ul style="list-style-type: none"> <li>• Albumin</li> <li>• Creatinine</li> <li>• Total Protein</li> <li>• Immunoglobulin A</li> <li>• Immunoglobulin M</li> </ul>	<u>Precinorm</u> <ul style="list-style-type: none"> <li>• Albumin</li> <li>• Creatinine</li> <li>• Total Protein</li> <li>• Immunoglobulin G</li> </ul> <u>Precipath</u> <ul style="list-style-type: none"> <li>• Albumin</li> <li>• Creatinine</li> <li>• Total Protein</li> <li>• Immunoglobulin A</li> <li>• Immunoglobulin M</li> <li>• Immunoglobulin G</li> </ul>

Note:

Page 31 is a duplicate of Page 30.

## 510(k) Summary, continued

Characteristic	Precinorm® PUC / Precipath® PUC (Predicate device, K040264)	Precinorm® PUC / Precipath® PUC (Modified Device)
Intended Use	For use in quality control by monitoring accuracy and precision for the quantitative methods as specified in the enclosed value sheet	Same
Format	Liquid ready-for-use control based on a buffered aqueous solution. Concentrations of control components have been adjusted to represent normal and pathological ranges.	Same
Stability	Unopened Stable at 2-8°C until expiration date Opened: Stable at 2 to 8°C for 4 weeks	Same
Constituent Analytes with Assigned Values	<u>Precinorm</u> <ul style="list-style-type: none"> <li>• Albumin</li> <li>• Creatinine</li> <li>• Total Protein</li> </ul> <u>Precipath</u> <ul style="list-style-type: none"> <li>• Albumin</li> <li>• Creatinine</li> <li>• Total Protein</li> <li>• Immunoglobulin A</li> <li>• Immunoglobulin M</li> </ul>	<u>Precinorm</u> <ul style="list-style-type: none"> <li>• Albumin</li> <li>• Creatinine</li> <li>• Total Protein</li> <li>• Immunoglobulin G</li> </ul> <u>Precipath</u> <ul style="list-style-type: none"> <li>• Albumin</li> <li>• Creatinine</li> <li>• Total Protein</li> <li>• Immunoglobulin A</li> <li>• Immunoglobulin M</li> <li>• Immunoglobulin G</li> </ul>



DEPARTMENT OF HEALTH & HUMAN SERVICES

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Theresa M. Ambrose, PhD, RAC  
Regulatory Affairs Principal  
Roche Diagnostics  
9115 Hague Road  
PO Box 50457  
Indianapolis, IN 46250

JAN 31 2005

Re: k050026  
Trade/Device Name: C.f.a.s (Calibrator for Automated Systems) Proteins in Urine/CSF (PUC)  
Precinorm® Proteins in Urine/CSF (PUC) and Precipath ® PUC  
Regulation Number: 21 CFR 862.1150  
Regulation Name: Calibrator  
Regulatory Class: Class II  
Product Code: JIX, JJY  
Dated: January 5, 2005  
Received: January 6, 2005

Dear Dr. Ambrose:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2 –

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific information about the application of labeling requirements to your device, or questions on the promotion and advertising of your device, please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (240)276-0484. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>

Sincerely yours,

A handwritten signature in cursive script that reads "Jean M. Cooper MS, DVM".

Jean M. Cooper, MS, D.V.M.

Director

Division of Chemistry and Toxicology

Office of *In Vitro* Diagnostic Device

Evaluation and Safety

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K050026

Device Name: Precinorm® Proteins in Urine/CSF (PUC) and Precipath® PUC

### Indications For Use:

Precinorm® PUC (Proteins in Urine/CSF) is for use in quality control by monitoring accuracy and precision for the quantitative methods as specified in the enclosed value sheet. Precipath® PUC (Proteins in Urine/CSF) is for use in quality control by monitoring accuracy and precision for the quantitative methods as specified in the enclosed value sheet.

Prescription Use X  
(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
Division Sign-Off

Office of In Vitro Diagnostic  
Device Evaluation and Safety

510(k) K050026

Page 1 of 1



## Indications for Use

510(k) Number (if known): K050026

Device Name: C.f.a.s. (Calibrator for Automated Systems) Proteins in Urine/CSF (PUC)

### Indications For Use:

C.f.a.s. PUC is for use in the calibration of quantitative Roche methods on Roche clinical chemistry analyzers as specified in the enclosed value sheet.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use             
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

  
Division Sign-Off

Office of In Vitro Diagnostic  
Device Evaluation and Safety

K050026

Page 1 of