

**SECTION VIII**

**510(K) SUMMARY OF SAFETY AND EFFECTIVENESS  
for 1cc Insulin SafePro\* Safety Syringe**

**1. REGULATORY AUTHORITY**

Safe Medical Device Act of 1990, CFR 807.92.

**2. CONTACT PERSON**

C. Y. Lin (Ms)  
Manager, Quality Assurance & Regulatory Affairs  
FORMOSA MEDICAL DEVICES, INC.

U.S. Liaison Office  
11497 Columbia Park Drive West, Suite #9  
Jacksonville, FL 32258

**3. NAME OF MEDICAL DEVICE**

Classification Name:	Syringe, Antistick
Classification Code:	MEG
Common/Usual Name:	Syringe
Proprietary Name:	1cc Insulin SafePro* Safety Syringe

**4. DEVICE CLASSIFICATION**

The General Hospital Panel has classified Antistick Syringes (21CFR880.5860) into Class II, Special Controls under section 513 of the Act.

**5. STATEMENT OF SUBSTANTIAL EQUIVALENCE**

The 1cc Insulin SafePro\* Safety Syringe is essentially equivalent to:

- A) B-D Conventional hypodermic and insulin syringe, K980580, and
- B) The 1cc SafePro Safety Syringe, K022063, in design, materials, instructions for use, and product claims.

**6. INTENDED USE**

The primary intended use for the 1cc Insulin SafePro\* Safety Syringe is aspirating and infusing insulin fluid.

The secondary intended use for the 1cc Insulin SafePro\* Safety Syringe is for needlestick protection; the device may aid in the reduction of needlestick injuries.

**7. DESCRIPTION OF DEVICE**

The 1cc Insulin SafePro\* Safety Syringe consists of a syringe assembly and a needle assembly. The device has a built-in safety feature to reduce the risk of accidental needlestick injuries. Other than scale graduation for use with U-100 strength insulin, the 1cc Insulin SafePro\* Safety Syringe is essentially identical to the legally marketed 1cc SafePro\* Safety Syringe in design and materials of construction.

**8. SUMMARY OF MATERIAL TESTING**

The materials of construction of the 1cc Insulin SafePro\* Safety Syringe are identical to those for the SafePro\* Safety Syringe and 1cc SafePro Safety Syringe. Those materials were already tested for material safety and biocompatibility as indicated in previous 510(K) submissions, K012726 and K022063. Therefore, no new biocompatibility tests are necessary.

**9. SUMMARY OF SIMULATED USE STUDY**

The safety feature and other functional and performance characteristics of the 1cc Insulin SafePro\* Safety Syringe are identical to those for the SafePro\* Safety Syringe and 1cc SafePro\* Safety Syringe. Those features and characteristics were already verified and validated as shown in the two previous 510(k) notifications K012726 and K022063. Therefore, no new tests are necessary.

**10. CONCLUSION**

The results obtained indicate that the 1cc Insulin SafePro\* Safety Syringe is safe and effective for its intended use.

### **COMPARISON MATRIX**

#### **1cc Insulin SafePro\* Safety Syringe, B-D Conventional Syringe, and 1cc SafePro\* Safety Syringe**

This matrix was developed in accordance with the ODE Guidance titled “Supplementary Guidance on the Content of Premarket Notification [510(K)] Submissions for Medical Device with Sharp Injury Prevention Features (March, 1995)”.

<b>Factor</b>	<b>1cc Insulin SafePro* Safety Syringe</b>	<b>B-D Conventional Insulin Syringe</b>	<b>1cc SafePro* Safety Syringe</b>
<b>510(K) Number</b>	To Be Assigned	K980580	K012726
<b>Intended Use and Claims</b>			
• Same Intended Use	Yes	Yes	Yes
• Tissue Puncture Device	Yes	Yes	Yes
• Sharps Injury Prevention Feature	Yes	No	Yes
• Reduces Risks of Accidental Needlesticks	Yes	No	Yes
• Delivers fluid to or draw blood from patient, short term	Yes	Yes	Yes
• Conventional Insertion Technique with minimal training	Yes	Yes	Yes
<b>Technological Features</b>			
• Same Technology Features	No	No	Yes
• Safety Mechanism - Manual Actuation	Yes	Not Applicable	Yes
• Safety Mechanism Remains Activated During Disposal	Yes	Not Applicable	Yes
• Safety Mechanism is an Integral Part of the Device	Yes	Not Applicable	Yes
<b>Material of Construction</b>			
• Needle • Needle Hub	• Stainless Steel •PP*	• Stainless Steel •PP	•Stainless Steel •PP* with more needle gauges (Different hub colorants)
• Barrel, Plunger, Sheath	PP	PP	PP
• Gasket	Synthetic Rubber	Synthetic Rubber	Synthetic Rubber

\*PP – Polypropylene resin



APR 5 2005

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. C.Y. Lin  
Manager, Quality Assurance & Regulatory Affairs  
FORMOSA MEDICAL DEVICES, INCORPORATED  
16F, No. 182, Section 2  
Tunhua South Road  
Taipei  
TAIWAN, R.O.C. 106

Re: K050134  
Trade/Device Name: 1cc Insulin\* SafePro Safety Syringe  
Regulation Number: 21 CFR 880.5860  
Regulation Name: Piston Syringe  
Regulatory Class: II  
Product Code: MEG  
Dated: January 21, 2005  
Received: January 24, 2005

Dear Ms. Lin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Chiu Lin', with a stylized flourish at the end.

Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known):

K050134

Device Name:

1cc Insulin SafePro\* Safety Syringe

Indications For Use:

The primary indicated use for the 1cc Insulin SafePro\* Safety Syringe is for aspirating and injecting insulin fluid.

The secondary indicated use for the 1cc Insulin SafePro\* Safety Syringe is for needlestick protection; the device may aid in the reduction of needlestick injuries.

Prescription Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use X  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER  
PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
Division Signatory  
Division of Anesthesiology, General Hospital,  
Injection Control, Dental Devices  
510(k) Number K050134

\*Trademark