

3.0 510(k) SummaryPage 1 of 1

Sponsor: Synthes (USA)
1302 Wrights Lane East
West Chester, PA 19380
(610) 719-5000

Device Name: Synthes 3.0mm Headless Compression Screws

Classification: 21 CFR 888.3040: Smooth or Threaded Metallic Bone Fixation Fastener.

Predicate Devices: Synthes 3.0 mm Sterile 3.0 mm Cannulated Screw and Threaded Washer

Device Description: The Synthes 3.0 mm Headless Compression Screws are cannulated and are self-drilling / self-tapping with a star-drive mechanism, and have a threaded head which can be countersunk into the bone. The screws are available in short and long thread lengths ranging from 10 mm to 40 mm. The screws are available in Stainless Steel and a Titanium Alloy.

Intended Use: The Synthes 3.0 mm Headless Compression Screws are intended for fixation of intra-articular and extra-articular fractures and non-unions of small bones and small bone fragments; arthrodeses of small joints; bunionectomies and osteotomies. Examples include, but are not limited to scaphoid and other carpal bones, metacarpals, tarsals, metatarsals, patella, ulnar styloid, capitellum, radial head and radial styloid.

Substantial Equivalence: Information presented supports substantial equivalence.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

APR 21 2005

Ms. Sheri L. Musgnung
Senior Regulatory Specialist
Synthes (USA)
1302 Wrights Lane
West Chester, Pennsylvania 19380

Re: K050636

Trade/Device Name: Synthes (USA) 3.0 mm Headless Compression Screws
Regulation Number: 21 CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: II
Product Code: HWC
Dated: March 10, 2005
Received: March 11, 2005

Dear Ms. Musgnung:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2- Ms. Sheri L. Musgnung

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0120 . Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Miriam C. Provost", is written over the typed name.

Miriam C. Provost, Ph.D.
Acting Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure



2.0

Indications for Use

510(k) Number (if known): _____

Device Name: Synthes (USA) 3.0 mm Headless Compression Screws

Indications for Use:

The Synthes 3.0 mm Headless Compression Screws are intended for fixation of intra-articular and extra-articular fractures and non-unions of small bones and small bone fragments; arthrodeses of small joints; bunionectomies and ostetotomies. Examples include, but are not limited to scaphoid and other carpal bones, metacarpals, tarsals, metatarsals, patella, ulnar styloid, capitellum, radial head and radial styloid.

Prescription Use X
(Per 21 CFR 801.109)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF
NEEDED)

(Division Sign-Off)

**Division of General, Restorative
and Neurological Devices**

Concurrence of CDRH, Office of Device Evaluation (ODE)

510(k) Number K050636