Astra Tech, Inc.

#### 510K Summary

# Astra Tech AB Lofric Primo Single Use Urinary Catheter

Sponsor Name
 Astra Tech, Inc.

#### 2. Device Name

**Proprietary Name:** 

Astra Tech AB LoFric® Primo™ Single Use

**Urinary Catheter** 

Common/Usual Name:

**Urethral Catheters** 

Classification Name:

Urethral Catheters and Accessories 72 GBM, Class II – Urology Devices

# 3. Identification of Legally Marketed Device

The modified LoFric Primo Single Use Urinary Catheter is substantially equivalent in intended use to the LoFric Single Use Urinary Catheter (K896750).

### 4. Device Description

The LoFric Primo Single Use Urinary Catheter is designed as an intermittent pathway for drainage of the bladder. The device consists of a catheter, coated with a hydrophilic low-friction coating.

The surface is hydrophilic and when the catheter is activated with the water integrated in the package, it becomes slippery and ready to use. The catheter is provided in a variety of lengths and sizes.

The catheter and the activation water is separated, sealed in one bag. By holding the product upright and exerting a light pressure on the folded water pocket, the water will run down and wet the catheter. The bag is opened and the catheter is inserted into the patient's urethra.[s2]

#### 5. Intended Use

The LoFric Primo Single Use Urinary Catheter is intended for intermittent catheterization of the urethra.

# 6 Comparison of Technological Characteristics

The modified The LoFric Primo Single Use Urinary Catheter is substantially equivalent in intended use and design to the currently marketed The LoFric Single Use Urinary Catheter (K896750).

The difference between the LoFric Primo Single Use Urinary Catheter and the predicate is the packaging concept, which includes the addition of water. The method of sterilization is changed. This difference does not raise new questions of safety and effectiveness. Laboratory data demonstrates this.

## 7. Performance Testing

Laboratory testing [s3] and biocompatibility testing was conducted to determine device functionality and conformance to design input requirements.



MAY - 5 2005

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Mr. Scott A. Root President Astra Tech, Inc. 890 Winter Street, Suite 310 WALTHAM MA 02451

Re: K050874

Trade/Device Name: LoFric® Primo™ Single Use Urinary Catheter

Regulation Number: 21 CFR §876.5130

Regulation Name: Urological catheter and accessories

Regulatory Class: II Product Code: GBM Dated: April 5, 2005 Received: April 7, 2005

Dear Mr. Root:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the <u>Code of Federal Regulations</u>, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

21 CFR 876.xxxx	(Gastroenterology/Renal/Urology)	240-276-0115
21 CFR 884.xxxx	(Obstetrics/Gynecology)	240-276-0115
21 CFR 892.xxxx	(Radiology)	240-276 <b>-</b> 0120
Other	(200001-85)	240-276-0100

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <a href="http://www.fda.gov/cdrh/dsma/dsmamain.html">http://www.fda.gov/cdrh/dsma/dsmamain.html</a>

Sincerely yours,

Nancy C. Brogdon

Nancy C. Brogdon

Director, Division of Reproductive, Abdominal, and Radiological Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

510(k) Number (if known):	K0508 14
Device Name: Astra	Tech AB LoFric® Primo™ Single Use Urinary Catheter
Indications For Use:	The LoFric® Primo™ Single Use Urinary Catheter is intended for intermittent catheterization of the urethra.
	<del></del>
(PLEASE DO NOT WRIT	TE BELOW THIS LINE. CONTINUE ON ANOTHER PAGE IF NEEDED)
Concurrence	e of CDRH, Office of Device Evaluation (ODE)
<i></i>	
Prescription Use(Per 21 CFR 801.109)	OR Over-The-Counter Use
	<i>A</i> • • •
D. ·	ision Sign-Off) sion of Reproductive, Abdominal,
510/L	Radiological Devices